 <p>Sirtex Medical Pty Ltd Shop 6, 207 Pacific Highway St Leonards NSW 2065 Australia Tel.: +61 2 9964 8400</p>	 <p>Sirtex Germany Manufacturing GmbH Branch Frankfurt am Main Industriepark Hoechst, building G808 65926 Frankfurt am Main Germany</p>	<p>Donawa Lifescience Consulting Srl Piazza Albania, 10 00153 Rome, Italy</p> <p>EC REP</p>	<p>MedEnvoy Switzerland Gotthardstrasse 28 6302 Zug Switzerland</p> <p>CH REP</p>	
--	--	---	--	---

SIR-Spheres® Y-90 resin microspheres

(Yttrium-90 Microspheres)

SIR-Y001

ENGLISH

1. DESCRIPTION

SIR-Spheres Y-90 resin microspheres consist of biocompatible microspheres containing yttrium-90. Yttrium-90 is a high-energy pure beta-emitting isotope with no primary gamma emission. The maximum energy of the beta particles is 2.27 MeV with a mean of 0.93 MeV. The half-life is 64.1 hours. The maximum range of emissions in tissue is 11 mm with a mean of 2.5 mm.

SIR-Spheres Y-90 resin microspheres are implanted into hepatic tumors by injection into the hepatic artery using a catheter. The SIR-Spheres Y-90 resin microspheres distribute non-uniformly in the liver, primarily due to the unique physiological characteristics of the hepatic arterial flow, the tumor to normal liver ratio of the tissue vascularity, and the size of the tumor. The tumor usually gets higher density per unit distribution of SIR-Spheres Y-90 resin microspheres than the normal liver. The density of SIR-Spheres Y-90 resin microspheres in the tumor is 5 to 6 times of the normal liver tissue. In general, 1 GBq of yttrium 90/kg of tissue provides 49.67 Gy of radiation dose. Once SIR-Spheres Y-90 resin microspheres are implanted into the liver, they are not metabolized or excreted, and they stay permanently in the liver with small amounts within the vasculature of normal parenchyma. The device is not phagocytized, nor does it dissolve or degrade after implantation. High dose radiation emitted from the device is cytotoxic to cells within the range of the radiation. SIR-Spheres Y-90 resin microspheres do not achieve primary intended purposes through chemical action within or on the body and do not depend upon being metabolized for the achievement of primary intended purposes and therefore have been designated a permanently implanted medical device. After the device has decayed, the non-radioactive microspheres remain intact and are not removed from the body.

The device has the potential to interact with other cytotoxic agents and is typically co-administered concomitantly with chemotherapeutic regimen.

This applies to chemotherapeutic agents used for the treatment of either the same tumors targeted by the microspheres, or distant metastases. This interaction may be exploited to the benefit of the patient, in that there can be an additive toxicity on tumor cells, which can enhance the tumor's cell kill rate.

As brachytherapy is defined as radiation treatment given by implanting radioactive devices directly into the body as close to the cancer as possible, SIR-Spheres Y-90 resin microspheres meet the definition of brachytherapy sources. Brachytherapy allows doctors to deliver higher doses of radiation to more-specific areas of the body, compared with the conventional form of radiation therapy (external beam radiation) that projects radiation from a machine outside of your body.

Each device is for single patient use.

2. INTENDED PURPOSE

SIR-Spheres Y-90 resin microspheres are intended for Brachytherapy and are permanently implanted via the hepatic artery to deliver therapy to hepatic tumors.

3. INDICATIONS FOR USE

SIR-Spheres Y-90 resin microspheres are indicated for the treatment of:

- unresectable hepatocellular carcinoma (HCC) or
- unresectable metastatic liver tumors from primary colorectal cancer (mCRC) in patients refractory to or intolerant of chemotherapy.

4. INTENDED PATIENT POPULATION

Patients with unresectable hepatocellular carcinoma (HCC) or unresectable metastatic liver tumors from primary colorectal cancer refractory to or intolerant of chemotherapy may be considered for treatment with SIR-Spheres. The ability to safely perform resection of tumor(s) is determined by review of individual patient characteristics, for instance during a multi-disciplinary tumor board (MDT). Patients with metastatic colorectal cancer (mCRC) are deemed refractory to or intolerant of chemotherapy after at least one drug or biologic. The potential benefit of selective internal radiation therapy for disease control within the liver is realized in patients with appropriate vascular anatomy of the liver and surrounding tissue, relatively good liver function or functional liver reserve, low lung shunt, and appropriate activity prescription.

5. INTENDED USERS

Intended users include Interventional Radiologists, Radiation Oncologists, and Nuclear Medicine Technicians trained in the use of Sirtex products. In addition, Radiation Safety Officers (RSO) are responsible to oversee appropriate use to ensure radiation safety policies are followed.

6. INTENDED USE ENVIRONMENT / CONDITIONS OF USE

The device is for professional use in the nuclear medicine lab setting and an angiography suite that are appropriately licensed to handle radioactive therapies.

7. MATERIALS OF CONSTRUCTION

Patients receiving SIR-Spheres Y-90 resin microspheres are in direct contact with the following materials:

- Yttrium 90 Chloride
- Cation Exchange Resin (Sodium Form)
- Water for Injection
- 99.99% Yttrium Oxide
- 0.1M Sulphuric Acid Solution
- Tribasic Phosphate
- Monobasic Phosphate

8. HOW SUPPLIED

SIR-Spheres Y-90 resin microspheres are provided in a vial with water for injection. Each vial contains 3 GBq of Y-90 ± 10% (at the time of calibration) in a total of 5 cc water for injection. Each vial contains

40 – 80 million microspheres with a diameter between 20 and 60 micrometers. The vial is shipped within a 6.4 mm minimum thickness lead pot. The package consists of a crimp sealed SIR-Spheres Y-90 resin microspheres glass vial within a lead pot, and a package insert within Type A package.

SIR-Spheres Y-90 resin microspheres are provided sterile using steam sterilization and are single use and shall not be re-sterilized.

Patient specific activities are prepared in accordance with Section 18.

The vial and its contents should be stored inside its transportation container at room temperature (15-25 °C; 59-77 °F) and a maximum relative humidity of 85%. Accessories should be stored in labeled product packaging at room temperature (15-25 °C; 59-77 °F) and a maximum relative humidity of 85%.

The calibration date (for radioactive contents) and the expiration information are printed on the vial label. The useful life of the SIR-Spheres Y-90 resin microspheres ends 24 hours after the time of calibration.

Inspect SIR-Spheres and accessories for any signs of damage or potential loss of sterility prior to use. SIR-Spheres Y-90 resin microspheres and accessories must not be used if the packaging is damaged or compromised. Replace with new and inspect SIR-Spheres Y-90 resin microspheres and accessories prior to procedure as needed. Report any compromised materials to Sirtex.

SIR-Spheres Y-90 resin microspheres + accessories (see Section 9) are single-use products (unless noted as a re-usable accessory) and must not be re-sterilized. Re-use and re-sterilization may create a risk of infection for the patient.

9. ACCESSORIES

SIR-Spheres Y-90 resin microspheres are intended to be used together with either Delivery System (accessories) listed below, for the dose preparation and implant procedure:

- Gen 1 Delivery System:
 - Gen 1 Delivery Box (SIR-B001)
 - Gen 1 Delivery Set (SIR-D001)
 - Gen 1 V-Vial (SIR-V001)
 - Gen 1 V-Vial Holder (SIR-H001)
 - Gen 1 Accessories Pack (SIR-K002)
 - Syringe Shield (SIR-S001)
- SIROS™ Delivery System:
 - SIROS Delivery Set (SIR-10300)
 - SIROS D-Vial Prep Set (SIR-10200)
 - SIROS Delivery Dome (SIR-10100)
 - Syringe Shield (SIR-S001)

The re-usable accessories (SIR-S001, SIR-B001 and SIR-10100) have been tested by Sirtex and shown to have a lifetime of a minimum of 100 functional cycles and 52 cleaning cycles as described per Section 24. Prior to use (and after cleaning per Section 24), the user must inspect the re-usable accessories for at least the following:

- Cracks, chips, or other visible damage to the components that negatively affects the use, safety or efficacy of the accessories
- Cracking, cloudiness, or discoloration of the components that negatively affects the visibility of the user
- Loss of functionality of components:
 - Inability to turn knobs
 - Lock/unlock features inoperable or damaged
 - Mating features damaged
 - Damage to threaded features

If any wear or damage is observed that may impact the function of the re-usable accessories, the user must discontinue use and contact Sirtex for replacements.

9.1 Components Common to both Delivery Systems

The Syringe Shield (SIR-S001) is intended to provide radiation shielding for the user(s) from SIR-Spheres (SIR-Y001) contained in the syringe during dose preparation. It is used with both the Gen 1 delivery system and the SIROS delivery system. The Syringe Shield is provided non-sterile and is a re-usable accessory (see Section 9 for details).

9.2 Gen 1 Delivery System

The Gen 1 Accessories Pack (SIR-K002) is intended to provide the Gen 1 Delivery Set (SIR-D001), Gen 1 V-Vial (SIR-V001) and Gen 1 V Vial Holder (SIR-H001). See individual accessories for details on sterility and use parameters.

The Gen 1 Delivery Box (SIR-B001) is intended to provide radiation shielding for the user(s) and patient, while allowing for the administration of SIR-Spheres (SIR-Y001). The Gen 1 Delivery Box is provided non-sterile and is a re-usable accessory (see Section 9 for details).

The Gen 1 V-Vial (SIR-V001) is intended to provide dose preparation accessories and to hold the prescribed dose of radioactive SIR-Spheres (SIR-Y001) for implantation. The Gen 1 V-Vial is provided in a sterile pouch using gamma sterilization and is a single use accessory and shall not be re-sterilized.

The Gen 1 V-Vial Holder (SIR-H001) is intended to provide radiation shielding for the user(s) and hold the V-Vial (SIR-V001) inside the Delivery Box (SIR-B001). The Gen 1 V-Vial holder is provided non-sterile and is a single use accessory and shall not be re-sterilized.

The Gen 1 Delivery Set (SIR-D001) is intended to be used inside the Delivery Box (SIR-B001) to allow the administration of SIR-Spheres (SIR Y001) from the V-Vial (SIR-V001) into the hepatic arterial circulation of the liver. The Gen 1 Delivery Set is provided in a sterile pouch using gamma sterilization and is a single use accessory and shall not be re-sterilized.

9.3 SIROS Delivery System

The SIROS Delivery Dome (SIR-10100) is intended to provide radiation shielding for the user(s) and patient, while allowing for the administration of SIR-Spheres (SIR-Y001). The SIROS Delivery Dome is provided non-sterile and is a re-usable accessory (see Section 9 for details).

The SIROS D-Vial Prep Set (SIR-10200) is held inside the Delivery Dome (SIR-10100) and is intended to provide dose preparation accessories, radiation shielding for the user(s), and to hold the prescribed dose of radioactive SIR-Spheres (SIR Y001) for implantation. The SIROS D-Vial Prep Set contains accessories provided non sterile and in sterile pouches using gamma sterilization and ethylene oxide sterilization (as applicable to the specific accessories). All accessories in the SIROS D-Vial Prep Set are single use and shall not be re-sterilized.

The SIROS Delivery Set (SIR-10300) is intended to be used inside the Delivery Dome (SIR-10100) to allow the administration of SIR-Spheres (SIR Y001) from the D-Vial (SIR-10200) into the hepatic arterial circulation of the liver. The SIROS Delivery Set is provided in a sterile pouch using gamma sterilization and is a single use accessory and shall not be re-sterilized.

Sections 19 - 24 provide additional safety precautions for users to ensure exposure to radiation is kept as low as reasonably achievable (ALARA).

10. CONTRAINDICATIONS

SIR-Spheres Y-90 resin microspheres are contraindicated in patients who have or are:

- had previous external beam radiation therapy to the liver
- ascites or are in clinical liver failure
- abnormal liver function tests (LFTs)
- total bilirubin > 2.0 mg/dL and/or albumin < 3.0 g/dL
- > 30 Gy absorbed radiation dose to the lungs, as estimated by the Technetium-99m macroaggregated-albumin (99mTc MAA) study described in Section 17
- pre-assessment angiogram that demonstrates potential pathways for deposition of microspheres to non-target organs such as the stomach, pancreas, or bowel that cannot be corrected (coiling)
- patients who are pregnant

11. WARNINGS

- Non-Target Delivery of SIR-Spheres Y-90 resin microspheres

Inadvertent delivery of SIR-Spheres Y-90 resin microspheres to extra-hepatic structures such as the stomach, duodenum, gallbladder, or pancreas may result in radiation injury to these structures, including but not limited to acute abdominal pain, acute gastritis, acute

cholecystitis, acute pancreatitis, and peptic ulceration. Angiographic technique must be employed to prevent the non-target delivery of SIR-Spheres Y-90 resin microspheres to any extra-hepatic structures.

• Radioembolization Induced Liver Disease (REILD)

Delivery of excessive radiation to the normal liver parenchyma may result in REILD. The risk of REILD may also be increased in patients with pre-existing liver disease. Consideration should be given to reducing the prescribed activity of SIR-Spheres Y-90 resin microspheres in the following clinical settings¹:

- Reduced liver functional reserve due to steatosis, steatohepatitis, hepatitis or cirrhosis
- Elevated baseline bilirubin level
- Prior hepatic resection
- Prior liver directed therapy
- Extensive prior treatment with systemic chemotherapy and/or biologic therapies

• Radiation Pneumonitis

High levels of implanted radiation and/or excessive shunting to the lung may lead to radiation pneumonitis. The lung radiation dose must be limited to ≤ 30 Gy for a single treatment session and ≤ 50 Gy cumulative dose.

• Other Tumor Types

The safety, performance, and benefit/risk profile of SIR-Spheres Y-90 resin microspheres in treating particular types of tumors outside the indications for use has yet to be established.

12. PRECAUTIONS

- Safety and effectiveness of this device in pregnant women, nursing mothers or children have not been established.
- A SPECT or PET scan of the upper abdomen be performed immediately after implantation of SIR-Spheres Y-90 resin microspheres. The SPECT or PET scan will detect radiation from the yttrium-90 to confirm placement of the microspheres in the liver.
- This product is radioactive. The use of this device is regulated, and local regulations must be followed when handling this device.
- Accepted radiation protection techniques should be used to protect staff when handling both the isotope and the patient.
- Patients may experience gastric problems following treatment, but proton pump inhibitors (PPI) or histamine H2-receptor antagonists (H2 blocking agents) may be used the day before implantation of SIR Spheres Y-90 resin microspheres and continued as needed to reduce gastric complications.
- Patients may experience abdominal pain immediately after administration of SIR Spheres Y-90 resin microspheres and pain relief may be required.
- SIR-Spheres Y-90 resin microspheres demonstrated a mild sensitization potential when tested dermally in an animal model.

13. ADVERSE EVENTS

When the patient is treated with proper technique, without excessive radiation to any organ, the common adverse events after receiving the SIR Spheres Y-90 resin microspheres are fever, transient decrease of haemoglobin, transient thrombocytopenia, mild to moderate abnormality of liver function tests (mild increase in aspartate aminotransferase, alkaline phosphatase, bilirubin), abdominal pain, nausea, vomiting, and diarrhea.

Potential Serious Adverse Events Due to High Radiation

Acute Pancreatitis: causes immediate severe abdominal pain. Verify by SPECT or PET of the abdomen and test for serum amylase.

Radiation Pneumonitis: causes excessive non-productive cough. Verify by X-ray evidence of pneumonitis.

Acute Gastritis: causes abdominal pain. Verify by standard methods to diagnose gastric ulceration.

Acute Cholecystitis: causes significant upper abdominal pain and may require cholecystectomy for resolution. Verify by appropriate imaging studies.

Radioembolization induced liver disease (REILD): REILD is a rare complication following Selective Internal Radiation Therapy (SIRT). REILD is characterized by a well-defined constellation of temporal, clinical, biochemical, and histopathologic findings. It typically manifests approximately

4 to 8 weeks post-SIRT and is characterized clinically by jaundice and ascites in the absence of tumor progression or bile duct obstruction.

The typical biochemical picture of REILD is an elevated bilirubin (> 3 mg/dL) in almost all cases, elevated alkaline phosphatase (ALP) and gamma- glutamyl transpeptidase (GGT) in most cases, accompanied by virtually no change in the transaminases (AST and ALT). If liver biopsy is performed, the typical histological appearance is of sinusoidal obstruction that may resemble veno-occlusive disease.

REILD may occur in both non-cirrhotic and cirrhotic patients.

Prophylactic treatment with methylprednisolone and ursodeoxycholic acid starting on the day of SIRT and continued for two months may reduce the incidence of REILD.

In the treatment of REILD, low molecular weight heparin may also be considered but both corticosteroids and heparin may only be useful if commenced very early in the course of the disease. See also Section 11 Warnings.

14. RESIDUAL RISK

Risks remaining after risk control measures have been implemented

- Unintended radiation exposure
- Delay to treatment
- Vascular system thrombosis/thrombus
- Infection
- Underdosing leading to suboptimal treatment

15. CLINICAL BENEFITS

The selective internal radiation therapy (SIRT) with SIR-Spheres Y-90 resin microspheres has the following clinical benefits:

- SIR-Spheres Y-90 resin microspheres efficiently treat HCC and mCRC liver tumors, measured in the form of improved (OS) and (DCR), over baseline.
- SIR-Spheres Y-90 resin microspheres safely deliver localized high dose radiation with optimal patient outcomes.
- SIR-Spheres Y-90 resin microspheres effectively and safely treat Hepatocellular Carcinoma (HCC) with vascular invasion (Portal Vein Tumor Thrombosis).

16. PERFORMANCE CHARACTERISTICS

The performance characteristics of SIR-Spheres Y-90 resin microspheres are related to:

- Appearance
- Radioactive content and properties
- pH levels
- Percentage of non-spherical particles
- Percentage of leachable Y-90
- Measurement of energy
- Single-use, steam sterilization
- Expiration of 24 hours post calibration
- Duration of contact with hepatic tissue
 - o SIR-Spheres Y-90 resin microspheres: > 30 days
 - o Y-90: Half-life of 64.1 hours

¹ Gil-Alzugaray et al. Hepatology, Vol 57, No. 3, 2013.

17. PRETREATMENT TESTING

17.1 Patient Tests Before Treatment with SIR-Spheres Y-90 resin microspheres

The following tests are performed before treatment to ensure the patient is suitable for treatment with SIR-Spheres Y-90 resin microspheres and determine the appropriate dose of SIR-Spheres Y 90 resin microspheres to prescribe:

- An hepatic angiogram to establish arterial anatomy of the liver
- ^{99m}Tc MAA Scan to determine the percent lung shunting to evaluate lung radiation dose exposure
- Biochemical tests of liver function
- CT, PET/CT, or MRI to determine the extent of the disease, lung, tumor and liver volumes for dosimetry calculations

17.2 Technique for Performing the Intra-Hepatic Technetium MAA Scan

To assess arterial perfusion of the liver and the fraction of radiopharmaceutical tracer that will pass through the liver and lodge in the lungs: Inject about 150 MBq of ^{99m}Tc MAA into the hepatic artery via a catheter.

Use a large field of view (FOV) gamma camera and obtain images of the thorax and abdomen (with the same acquisition time).

Draw Region of Interest (ROI) around the whole liver and the whole lung and get the total counts for the lung and the liver.

17.3 Lung Shunt Calculation Procedure

Calculate the lung shunt fraction (L) using following formula:

$$L = \left(\frac{\text{Lung Counts}}{\text{Liver Counts} + \text{Lung Counts}} \right) \quad \text{Equation 1}$$

In order to optimize the risk versus benefit for patients receiving SIR-Spheres Y-90 resin microspheres, limiting the radiation exposure to the lungs to < 30 Gy is required. The calculation of estimated radiation exposure to the lungs is given by the following formula:

Activity that may potentially reach the lung:

$$A_{\text{Lung}} = A_{\text{total}} \times L \quad \text{Equation 2}$$

Where:

A_{Lung} = lung activity [GBq]

A_{total} = total prescribed activity [GBq]

L = lung shunt fraction

The resulting lung dose, given that a given amount of activity shunts from the liver to the lung:

$$D_{\text{Lung}} = \frac{49670 \times A_{\text{Lung}}}{M_{\text{Lung}}} \quad \text{Equation 3}$$

Where:

D_{Lung} = lung dose [Gy]

A_{Lung} = lung activity [GBq]

M_{Lung} = mass of the lung [g]

18. CALCULATION OF INDIVIDUAL DOSE

There are two accepted methods for calculating the patient radiation dose, these being the Body Surface Area (BSA) model and the partition model.

18.1 The Body Surface Area (BSA) Method

The BSA method varies yttrium-90 activity according to the size of the patient and the size of the tumor within the liver. The BSA method can be used when treating the volume of a single treated lobe, as well as that of the entire liver. The approach of lobar treatment versus whole liver treatment with SIR-Spheres Y-90 resin microspheres is based on the presence of visible

tumors on pre-treatment CT or MR imaging. If liver tumors are only visible in one lobe, then SIR Spheres Y-90 resin microspheres should be administered to that lobe only, thus sparing the contralateral lobe from unnecessary internal radiation.

BSA must first be determined and is calculated from the following equation:

$$BSA = 0.20247 \times H^{0.725} \times W^{0.425} \quad \text{Equation 4}$$

Where:

BSA = Body Surface Area

H = height in meters

W = weight in kilograms

18.1.1. BSA prescribed activity calculation for whole liver / bilobar treatment

$$A_{\text{Admin}} = (BSA \cdot 0.2) + \left(\frac{V_{\text{Tumor}}}{V_{\text{Tumor}} + V_{\text{normal liver}}} \right) \quad \text{Equation 5}$$

Where:

A_{Admin} = SIR-Spheres activity to implant [GBq]

V_{Tumor} = volume of tumor

$V_{\text{normal liver}}$ = volume of non-tumor tissue in the treated volume

18.1.2. BSA prescribed activity calculation for lobar or super-selective treatment

In patients who receive lobar or segmental treatment with SIR-Spheres Y-90 resin microspheres, the prescribed activity must be reduced in accordance with the size of the portion of the liver being treated.

$$A_{\text{Admin}} = \left[BSA - 0.2 + \left(\frac{V_{\text{Tumor}}}{V_{\text{Treated}}} \right) \right] \times \left[\frac{V_{\text{Treated}}}{V_{\text{Liver}}} \right] \quad \text{Equation 6}$$

Where:

A

$Admin$ = SIR-Spheres activity to implant [GBq]

V_{Tumor} = volume of tumor in the portion of the liver being treated (i.e., lobe)

V_{Treated} = volume of the portion of the liver being treated AND the tumor volume

V_{Liver} = total volume of the whole liver, including tumor

BSA = Body Surface Area as per Equation 4

18.2 Partition Model

This method involves selecting safe radiation doses to the normal liver and lung and implanting the maximum activity that will not exceed these limits. The dose received by the tumor has no upper limit.

The partition model should be used where the tumor mass is a discrete area within the liver. The technique requires two measurements to be made:

1. Measurement of the volume of tumor and normal perfused liver determined from a CT or MR scan.
2. Measurement of the proportions of ^{99m}Tc MAA activity that lodges in the tumor, normal perfused liver and lung as determined from imaging.

Therefore, to calculate the activity to be implanted, it is necessary to:

- Determine the volumes of the normal perfused liver and tumor and convert each volume to mass²
- Determine the volume of the lung and convert to mass³
- Using the ^{99m}Tc MAA scan, determine the activity in the lung, tumor and perfused normal liver

² ICRP, 2015. Radiation Dose to Patients from

Radiopharmaceuticals: A Compendium of Current Information Related to Frequently Used Substances. ICRP Publication 128. Ann. ICRP 44(2S).

³ Kappadath SC, Lopez BP, Salem R, Lam MG. Lung shunt and lung dose calculation methods for radioembolization treatment planning. Q J Nucl Med Mol Imaging.

- Determine the Tumor to Normal Activity Ratio (TNR) calculated as activity per unit mass of organ or tissue using the following equation:

$$TNR = \frac{A_{\text{Tumor}}/M_{\text{Tumor}}}{A_{\text{PNL}}/M_{\text{PNL}}} = \frac{\text{average counts or activity in tumor}}{\text{average counts or activity in perfused normal liver}} \quad \text{Equation 7}$$

Where:

A_{Tumor} = Activity in tumor

M_{Tumor} = Mass of tumor

A_{PNL} = Activity in the perfused normal liver

M_{PNL} = Mass of the perfused normal liver

To calculate the total activity to be implanted, use the following equations. The activity required should be calculated to accommodate the lung and normal tissue doses as the limiting factors.

$$A_{\text{Admin}} = \frac{D_{\text{PNL}} (TNR * M_{\text{Tumor}} + M_{\text{PNL}})}{49670 * (1-L)} \quad \text{Equation 8}$$

Where:

A_{Admin} = SIR-Spheres activity to implant [GBq]

D_{PNL} = Mean absorbed dose to perfused normal liver [Gy]

M_{Tumor} = Mass of tumor [g]

M_{PNL} = Mass of perfused normal liver [g]

L = Lung Shunt Fract

and,

$$D_{\text{Tumor}} = TNR * D_{\text{PNL}} \quad \text{Equation 9}$$

Where:

D_{Tumor} = Desired absorbed dose to tumor [Gy]

19. GEN 1 DOSE PREPARATION PROCEDURE (IFU-VV current revision)

19.1 Gen 1 Accessories provided by Sirtex:

- Gen 1 V-Vial (SIR-V001)
- Gen 1 V-Vial Holder (SIR-H001)
- Syringe Shield (SIR-S001)

19.2 Additional supplies needed:

- Sterile water for Injection
- Two 25 G vent needles with filters
- 5 mL syringe
- Alcohol swabs
- Forceps
- Dose Calibrator (ion chamber)
- Radiation survey meter
- Radiation spill kit

19.3 Gen 1 Dose Prep Procedure

- a. Unpack SIR-Spheres Y-90 resin microspheres, leaving the glass shipping vial in its lead pot and place on the benchtop.
- b. Remove the Gen 1 V-Vial from its sterile pouch, then remove the center of the aluminum seal with forceps and wipe the rubber septum with an alcohol swab.
- c. Place the Gen 1 V-Vial in the Gen 1 V-Vial Holder and screw the Gen 1 V-Vial Holder Lid for stability and shielding.
- d. Insert a short 25 G needle with filter through the rubber septum of the Gen 1 V-Vial until it just pierces the septum to create a vent. It is recommended that either a purpose-designed venting needle or a separate filter attached to the short 25 G needle is used in order to prevent fluid leakage.
- e. Leaving the SIR-Spheres Y-90 resin microspheres glass shipping vial in the lead pot, invert the lead pot and shake vigorously before opening to re-suspend the SIR-Spheres Y-90 resin microspheres, which will have settled during shipping. Resuspension ensures a homogenous solution for dose preparation.
- f. Quickly open the lead pot and remove the glass shipping vial using forceps.
- g. Determine the total activity of SIR-Spheres

Y-90 resin microspheres in the glass shipping vial using an appropriate ion chamber (dose calibrator) and then return the glass shipping vial to the lead pot.

- h. Determine the volume of the SIR-Spheres Y-90 resin microspheres that needs to be withdrawn from the glass shipping vial to provide the intended patient specific activity of SIR-Spheres Y-90 resin microspheres.
 - i. Partially remove the center of the aluminium seal of the SIR-Spheres Y-90 resin microspheres glass shipping vial and wipe the rubber septum with an alcohol swab.
 - j. Insert a 25 G needle with filter through the rubber septum of the glass shipping vial to create a vent, ensuring the needle is well clear of the contents in the glass shipping vial. It is recommended that either a purpose-designed venting needle or a separate filter attached to the short 25 G needle is used in order to prevent fluid leakage.
 - k. Attach a lubricated 21 G needle at least 50mm in length to a 5 mL Luer lock syringe and place it into the Syringe Shield.
 - Unscrew the top of the Syringe Shield and place the 5 mL syringe inside.
 - Place the top of the Syringe Shield back on and tighten securely.
 - l. Using the Syringe Shield and lubricated 21 G needle, puncture the rubber septum of the SIR-Spheres Y-90 resin microspheres glass shipping vial and quickly draw back and forth at least 6 times in order to re suspend the SIR-Spheres Y-90 resin microspheres thoroughly.
 - m. Quickly withdraw the volume of SIR Spheres Y-90 resin microspheres that will provide the intended patient specific activity, being careful of needlestick injuries, re-cap the needle safely and set aside Syringe Shield and needle in prepared radiation work area
 - n. Gently agitate the glass shipping vial to re dispense the microspheres and measure the activity remaining in the glass shipping vial with the dose calibrator. Subtract the activity remaining in the glass shipping vial from the starting total activity in the glass shipping vial to determine the amount of activity that has been drawn up into the 5 mL syringe.
 - o. If the amount of activity that has been drawn up into the 5 mL syringe is not correct, transfer the SIR-Spheres Y-90 resin microspheres back into the glass shipping vial and redraw the necessary volume of SIR-Sphere Y-90 resin microspheres.
- Once the correct activity has been drawn up, transfer the SIR-Spheres Y-90 resin microspheres from the 5 mL syringe into the vented Gen 1 V-Vial in the Gen 1 V-Vial Holder. If the total volume in the 5 mL syringe is less than 3 mL, draw up enough sterile water for injection to make up to a total volume of 3 – 5 mL before transferring the SIR-Spheres Y-90 resin microspheres into the V-Vial. Ensure that the distance between any puncture holes in the rubber septum of the Gen 1 V-Vial are at least 2mm apart. CAUTION: This step should be done ONCE only.
- p. Remove the vent needle from the Gen 1 V Vial. Ensure the screw on lid of the Gen 1 V-Vial Holder is secure and place (do not force) the black plug into place.
 - q. Remove the vent needle from the glass shipping vial and replace the lid on the lead pot.
 - r. The patient specific activity of SIR-Spheres Y-90 resin microspheres is now ready for transport to the angiography suite in which the implantation will be performed.

20. GEN 1 IMPLANT PROCEDURE (IFU-GLO current revision)

Doctors must refer to the Sirtex Medical Pty Ltd

Training Manual (TRN-EU-001) for delivering SIR Spheres Y-90 resin microspheres before attempting to implant this device.

20.1 Gen 1 Accessories provided by Sirtex:

- Gen 1 Delivery Set (SIR-D001)
- Gen 1 Delivery Box (SIR-B001)
- Gen 1 V-Vial (SIR-V001)
- Gen 1 V-Vial Holder (SIR-H001)

20.2 Additional supplies needed:

- Two 20 mL Luer Lock Syringes filled with non ionic solution (either 5% glucose/dextrose or water for injection).
- Do not use saline
- Patients with diabetes should be monitored for hyperglycemia if 5% dextrose/glucose is used to implant SIR Spheres Y-90 resin microspheres
- One 20 mL Luer Lock Syringe filled with non- ionic contrast
- A microcatheter with an inner diameter of at least 0.021" and a 45° tip configuration for vessels with a 90° angle of origin.
- Sterile absorbent material for protection of angiography suite
- Alcohol swabs
- Forceps or hemostat
- Radiation survey meter
- Radiation spill kit

20.3 GEN 1 PROCEDURE

The hepatic artery catheter is inserted by a trained interventional radiologist under x-ray guidance. This method allows complete control of exactly where the catheter is placed and allows routine checking of the catheter position throughout the implant procedure.

It is essential that SIR-Spheres Y-90 resin microspheres are not delivered to other organs, in particular the pancreas, stomach, or duodenum. If there is any possibility of SIR-Spheres Y-90 resin microspheres passing down the Gastroduodenal artery (GDA) then the implantation must not proceed. It may be preferable to block shunting vessels with an intraluminal coil or other agent to prevent SIR-Spheres Y-90 resin microspheres from flowing to non-target organs.

Note: Virtually all complications from SIR-Spheres Y-90 resin microspheres arise from the inadvertent delivery of SIR-Spheres Y-90 resin microspheres into small blood vessels that go to the pancreas, stomach or duodenum.

The radiologist must repeatedly check the position of the catheter during the procedure to ensure it remains correctly sited and that reflux of the SIR Spheres Y-90 resin microspheres into other organs does not occur. This is performed by injecting contrast medium through the B line port of the Gen 1 Delivery Set during the delivery of SIR-Spheres Y-90 resin microspheres. Contrast must not be administered on the D line port.

SIR-Spheres Y-90 resin microspheres must be delivered slowly at a rate of no more than 5 mL per minute. Rapid delivery may result in a more concentrated suspension of SIR-Spheres Y-90 resin microspheres that can cause clogging in the microcatheter or reflux back down the hepatic artery and into other organs. At the conclusion of the procedure, the catheter is removed.

- a) Firmly place Gen 1 Delivery Set 3-way stopcock into bracket on back wall of Gen 1 Delivery Box.
- b) Holes in the Gen 1 Delivery Box are color coded and marked with 'A', 'B' and 'D'. From inside the Gen 1 Delivery Box, insert Gen 1 Delivery Set lines through the corresponding holes.
- c) Remove caps from the end of the Gen 1 Delivery Set 'B' and 'D' lines and attach 20 mL syringes with a non-ionic solution (NOT SALINE). Prime all lines with non-ionic solution. This is done with covers left on needles to maintain sterility.
- d) In order to prime the 'C' line of the Gen 1 Delivery Set, while still engaged with the off handle, pull the Gen 1 Delivery Box

Stopcock Control Knob just past the limiting notch and rotate Control Knob 90 degrees (one-quarter turn) counterclockwise (left).

- Non-ionic solution will then flow through the 3-way stopcock into the 'C' line of the Gen 1 Delivery Set.
 - Re-engage the Stopcock Control Knob fully so it is limited to one-quarter turn.
- e) Place Gen 1 V-Vial Holder containing Gen 1 V Vial with SIR-Spheres Y-90 resin microspheres into the Retaining Ring in the Gen 1 Delivery Box. Remove the black plug from the Gen 1 V-Vial Holder. Swap Gen 1 V Vial septum with alcohol, e.g., skin wipes. Care must be taken when inserting needles to not contaminate them. If contamination occurs, then discard the Gen 1 Delivery Set and obtain a new one.
 - f) Insert Gen 1 Delivery Set 'C' Line needle through center of the Gen 1 V-Vial septum. Take care not to scrape the sidewalls of the Gen 1 V-Vial.
 - g) Insert Gen 1 Delivery Set 'D' Line needle through one side of the septum until needle rest at base of the 'V'. Ensure spacing between any punctured holes in the Gen 1 V-Vial septum is at least 2 mm. It is important that the Gen 1 Delivery Set 'D' Line needle goes to the bottom of the Gen 1 V-Vial to ensure that SIR Spheres Y-90 resin microspheres will swirl into suspension when non-ionic solution is injected. The SIR-Spheres Y-90 resin microspheres will then be decanted from the top of the solution. An excessively concentrated suspension of SIR-Spheres Y 90 resin microspheres may cause clogging in the microcatheter.
 - h) Remove cap on Gen 1 Delivery Set 'A' Line and connect to patient's microcatheter. After connecting, contrast media can be administered via Gen 1 Delivery Set 'B' Line by turning the Gen 1 Delivery Box Stopcock Control Knob to "Flush/Contrast" position. DO NOT USE ionic contrast media with SIR Spheres Y-90 resin microspheres directly.
 - i) When the Gen 1 Delivery System is fully assembled and the lid of the Gen 1 Delivery Box is in place, injecting non-ionic solution from 20 mL syringe on the Gen 1 Delivery Set 'D' Line will cause SIR-Spheres Y-90 resin microspheres to re-suspend. Turn the Gen 1 Delivery Box Stopcock Control Knob to "SIR Spheres" position and slowly deliver SIR- Spheres Y-90 resin microspheres at a rate of approximately 5 mL per minute. In order to achieve a slow and controlled delivery rate and maintain SIR-Spheres Y-90 resin microspheres in suspension, flow rate from the syringe may be given in pulses of 0.25 mL –0.5 mL. Use all 20 mL of non-ionic solution.
 - j) Avoid prolonged stoppages or interruptions to administration process, as this may allow SIR Spheres Y-90 resin microspheres to settle out of suspension, making them difficult to resuspend and potentially cause blockages.
 - k) Pulse the Gen 1 Delivery Set 'D' line syringe to slowly bring SIR-Spheres Y-90 resin microspheres into dilute suspension.
 - l) Avoid sideways tension Gen 1 Delivery Set tubing to avoid potential leaks at needle entry to rubber septum of Gen 1 V-Vial.
 - m) Do not over-pressurize the system, such as could occur if using syringes smaller than those recommended.
 - n) There will be some solution and SIR-Spheres Y-90 resin microspheres left in the Gen 1 V Vial once all non-ionic solution in the 20 mL syringe connected to the Gen 1 Delivery Set 'D' Line has been delivered. In order to deliver the remaining activity, without removing the Gen 1 Delivery Set 'D' Line needle, raise it up 15-20 mm and carefully push the Gen 1 Delivery Set 'C' Line needle to the bottom of the Gen 1 V-Vial. The remaining solution will empty

from the Gen 1 V-Vial by injecting air into the 'D' Line (Approx. 8-10 mLs). Care must be taken to prevent air from entering tubing going to the patient.

- o) Upon completion, follow the disposal and cleaning procedures in Section 24.

21. SIROS DOSE PREPARATION PROCEDURE (IFU-002-US current revision)

21.1 SIROS Accessories provided by Sirtex:

- SIROS D-Vial Prep Set (SIR-10200)
- Syringe Shield (SIR-S001)

21.2 Additional supplies needed:

- Sterile water for Injection
- 5 mL syringe
- 20 mL syringe
- Alcohol swabs
- Forceps
- Dose Calibrator (ion chamber)
- Radiation survey meter
- Radiation spill kit

21.3 SIROS Unpacking Procedure

- a) Unpack SIR-Spheres Y-90 resin microspheres lead pot from the plastic shipping pail and place it on the benchtop.
- b) Insert the SIROS D-Vial Holder into the SIROS Transport Base and place it on the prep surface area.
- c) Remove SIROS D-Vial from sterile pouch and place in SIROS D-Vial Holder.
- d) Remove the aluminum seal from the center of the SIROS D-Vial with forceps and wipe the rubber septum with an alcohol swab.

21.4 SIROS D-Vial Priming Procedure

- e) Connect the 0.22 micron filter to the 25 G vent needle. Insert it through the rubber septum of the SIROS D-Vial to create a vent.
- f) Ensure that vent needle tip is above the fill level marked on the SIROS D-Vial.
- g) Remove the blue cap from the SIROS D-Vial 'D' line.
- h) Attach a syringe filled with at least 10 mL of D5W or sterile water for injection to prime the SIROS D-Vial 'D' line. Fill the SIROS D-Vial to fill level marking. Ensure there is no air in the 'D' line.
- i) Disconnect priming syringe from the SIROS D-Vial 'D' line. Attach a new blue cap to the 'D' line connector.
- j) Remove blue cap from the SIROS D-Vial 'C' line.
- k) Connect same priming syringe to SIROS D-Vial 'C' line.
- l) Slowly draw fluid into the syringe from the SIROS D-Vial to prime the 'C' line until fluid is at the prime level marking on the SIROS D-Vial. Ensure there is no air in the 'C' line. CAUTION: Do not draw below the prime level marking.
- m) Disconnect syringe from the SIROS D-Vial 'C' line. Attach a new blue cap to the 'C' line connector.

21.5 SIROS Dose Draw Procedure

- n) Invert the lead pot and shake vigorously before opening to re-suspend the SIR Spheres Y-90 resin microspheres, which will have settled during shipping.
- o) Quickly open the lead pot and remove the shipping vial using forceps.
- p) Determine the total activity of SIR-Spheres Y-90 resin microspheres in the shipping vial using an appropriate ion chamber (dose calibrator) and then return the shipping vial to the lead pot.
- q) Determine the volume of SIR-Spheres Y-90 resin microspheres suspension that needs to be withdrawn from the glass shipping vial to provide the intended patient-specific activity of SIR-Spheres Y 90 resin microspheres.
- r) Partially peel back the aluminum seal of the SIR-Spheres Y-90 resin microspheres glass shipping vial and wipe the septum with an alcohol swab.
- s) Insert a 25-gauge needle through the

septum of the shipping vial to create a vent, ensuring the needle is well clear of the contents in the shipping vial.

- t) Connect the 0.22 micron filter to the 25 G vent needle. Insert through the rubber septum of the glass shipping vial to create a vent. Ensure that the tip of the needle is well clear of the contents.
- u) Attach the 21 G drawing up needle to a 5 mL Luer lock syringe and place it in the acrylic Syringe Shield.
- v) Use the shielded syringe and 21 G needle, puncture the septum of the SIR-Spheres Y 90 resin microspheres glass shipping vial, and quickly draw back and forth at least 6 times in order to re-suspend the SIR- Spheres Y-90 resin microspheres thoroughly.
- w) Quickly withdraw the determined volume of SIR-Spheres Y-90 resin microspheres suspension that will provide the intended patient specific activity. Prior to removing the syringe needle from the glass shipping vial, draw some air through the needle into syringe to draw spheres within needle up into the shielded syringe.
- x) Withdraw the needle from the glass shipping vial septum and re-cap the needle using forceps. Set aside on the benchtop.
- y) Using forceps, swirl the glass shipping vial to re-suspend the microspheres and measure the activity remaining in the glass shipping vial with the dose calibrator.
- z) Replace the glass shipping vial into the lead pot.
- aa) Subtract activity remaining in the glass shipping vial from the starting total activity to determine the amount of activity that has been drawn up into the 5 mL syringe.
- bb) If the amount of the activity that has been drawn up into the 5mL syringe is not correct, transfer the SIR-Spheres Y-90 resin microspheres suspension back into the glass shipping vial and repeat the previous steps to obtain the prescribed level of activity.
- cc) Once the correct activity has been obtained, remove the vent needle from the glass shipping vial and secure the lead pot cover.
- dd) Transfer the SIR-Spheres Y-90 resin microspheres from the 5mL syringe into the vented SIROS D-Vial. CAUTION: Do not exceed the fill level mark on the SIROS D-Vial.
- ee) Remove all needles from the SIROS D-Vial septum.
- ff) Screw on the SIROS D-Vial Holder cap.
- gg) The patient specific activity of SIR-Spheres Y-90 resin microspheres is now ready for transport to the angiography suite in which the implantation procedure will be performed.

22. SIROS IMPLANT PROCEDURE (IFU-001-US current revision)

Doctors must refer to the Sirtex Medical Pty Ltd Training Manual (TRN-EU-001) for delivering SIR Spheres Y-90 resin microspheres before attempting to implant this device.

22.1 SIROS Accessories provided by Sirtex:

- SIROS Delivery Dome (SIR-10100)
- SIROS D-Vial Prep Set (SIR-10200)
- SIROS Delivery Set (SIR-10300)

22.2 Additional supplies needed:

- Two 20 mL Luer Lock Syringes filled with non-ionic solution (either 5% glucose/dextrose or water for injection):
 - Do not use saline
 - Patients with diabetes should be monitored for hyperglycemia if 5% dextrose/glucose is used to implant SIR Spheres Y-90 resin microspheres
- One 20 mL Luer Lock Syringe filled with non- ionic contrast
- Microcatheter: the microcatheters should have an inner diameter of at least 0.021"

and a 45° tip configuration for vessels with a 90° angle of origin.

- Sterile absorbent material for protection of angiography suite
- Alcohol swabs
- Forceps or hemostat
- Radiation survey meter
- Radiation spill kit

22.3 SIROS Implant Procedure Setup

- a) Place the SIROS Delivery Dome on a stable table or cart near the patient. Visually inspect SIROS before use, do not use if damage or wear that could impact function is observed. Ensure that the end marked 'A' is directed toward the microcatheter.
- b) Remove the lid from the base by pressing the black unlock button. Set aside lid.

22.4 SIROS Delivery Set Setup

- c) Remove SIROS Delivery Set from sterile pouch.
- d) Position SIROS Delivery Set on the SIROS Delivery Dome with the 'A' line directed toward the microcatheter.
- e) Secure SIROS Delivery Set lines into their corresponding cleats on the SIROS Delivery Dome.

22.5 SIROS Priming Lines:

- f) Fill two 20 mL syringes with non-ionic solution.
- g) Remove the cap from the SIROS Delivery Set 'D' line.
- h) Connect on syringe to the SIROS Delivery Set 'D' line and inject non-ionic solution until all air is removed from the 'D' line and fluid exits the 'A' line.
- i) Fill two 20 mL syringes with non-ionic contrast solution.
- j) Remove the cap from the SIROS Delivery Set 'B' line.
- k) Connect one syringe to the SIROS Delivery Set 'B' line and inject until all air is removed for the 'B' line and exits the 'A' line.
- l) Remove vented cap from the SIROS Delivery Set 'A' line and attach the 'A' line to the microcatheter that is in position in the patient.
- m) CAUTION: Ensure the patient microcatheter is fully connected to the SIROS Delivery Set 'A' line before proceeding to step 22.6.

22.6 SIROS D-Vial Connection

- n) Remove the SIROS D-Vial Holder from its Transport Base. Place the D-Vial Holder into the cupped area on the lower shelf of the SIROS Delivery Dome and ensure it is fully seated.
- o) Remove the cap from the SIROS D-Vial 'D' line.
- p) Disconnect the SIROS Delivery Set 'D' line from the red connector and connect it to the SIROS D-Vial 'D' line.
- q) Remove the cap from the SIROS D-Vial 'C' line.
- r) Remove and discard the red connector from the SIROS Delivery Set 'C' line. Ensure a wet-to-wet connection is maintained on the 'C' line fittings.
- s) Connect the male Luer fitting on the SIROS Delivery Set 'C' line to the SIROS D-Vial 'C' line.
- t) Place the lid back onto the SIROS Delivery Dome base and press down to lock in place. Ensure SIROS Delivery Set tubing lines are not pinched

22.7 SIROS Delivery Procedure

- u) Deliver the SIR-Spheres Y-90 resin microspheres in approximately 2 mL aliquots at approximately 5 mL per minute infusion rate using the SIROS Delivery Set 'D' line syringe.
- v) Using the SIROS Delivery Set 'B' line syringe, follow each 2 mL aliquot of SIR Spheres Y-90 resin microspheres

suspension with a 2 mL aliquot of non-ionic contrast media to facilitate fluoroscopic visualization of flow.

- w) The delivery has been completed when at least 20 mL of non-ionic solution has been administered through the SIROS D-Vial.
- x) Upon completion, follow the disposal and cleaning procedures in Section 24.

23. DISASSEMBLY

23.1 Gen 1 Disassembly

- a) Once the infusion is complete, remove lid from the Gen 1 Delivery Box.
- b) Do not disconnect the patient catheter from the Gen 1 Delivery Set.
- c) Carefully remove catheter from the patient, while still attached to the Gen 1 Delivery Set. Use care when removing and handling the catheter as it may be contaminated with radioactivity. Coil the catheter and wrap within a sterile towel.
- d) Leaving everything attached, use a hemostat to push down the needles into the shielded Gen 1 V-Vial and place it (along with all potentially contaminated procedural accessories) into the designated receptacle for radioactive material and for post implant measurements.
- e) Remove 5 mL syringe from Syringe Shield and place it into the designated receptacle for radioactive material. Use care when removing and handling the syringe as it may be contaminated with radioactivity and for post implant measurements.
- f) Radioactive materials must be stored in accordance with local regulations governing radioactive materials storage.
- g) Ensure that all items that have come in contact with SIR-Spheres Y-90 resin microspheres are appropriately disposed of as they may be radioactive.

23.2 SIROS Disassembly

- a) Once the infusion is complete, remove syringes from SIROS Delivery Set 'B' and 'D' lines. Do not disconnect the patient catheter from the SIROS Delivery Set.
- b) Remove the lid from the base of the SIROS Delivery Dome by pressing the black unlock button. Set aside lid on the stable table or cart.
- c) Carefully remove catheter from the patient, while still attached to the SIROS Delivery Set. Use care when removing and handling the catheter as it may be contaminated with radioactivity. Coil the catheter and wrap with a sterile towel.
- d) Leaving everything attached, lift the SIROS D-Vial Holder out of the cupped area of the SIROS Delivery Dome and place it (along with all potentially contaminated procedural accessories) into the waste container for radioactive material and for post implant measurements.
- e) Remove 5 mL syringe from Syringe Shield and place it into the designated receptacle for radioactive material. Use care when removing and handling the syringe as it may be contaminated with radioactivity and for post implant measurements.
- f) Radioactive materials must be stored in accordance with local regulations governing radioactive materials storage.
- g) Ensure that all items that have come in contact with SIR-Spheres Y-90 resin microspheres are appropriately disposed of as they may be radioactive.

24. CLEANUP, WASTE DISPOSAL AND STORAGE

- Following dose preparation and delivery, survey equipment and people using appropriate radiation survey meter to ensure any contamination is contained.
- If re-usable accessories (i.e., Syringe Shield or Gen 1 Delivery box) are visibly soiled, clean surfaces by spraying with cleaning

solution (0.5% sodium hypochlorite, NaOCl) until thoroughly wet. Allow it to stand for 1 minute. Wipe with a soft cloth until all visible soil is removed.

- To disinfect, spray all surfaces with 0.5% NaOCl solution at 6-8 inches until thoroughly wet. Allow to stand for at least 1 minute. Then wipe with a soft cloth to remove any residue.
- For cleaning and low-level disinfection of SIROS Delivery Dome, use 4 each of a quaternary / alcohol based cleaning / disinfecting wipe such as Super Sani-Cloth. Wipe entire exposed surface and visually inspect to ensure surfaces are clean. Allow surfaces to remain wet for duration described in the cleaning / disinfecting wipe instructions. Dry surfaces with a soft cloth as needed. Do NOT spray liquid cleaners or use abrasive materials on SIROS Delivery Dome.
- Do not use alcohol-based cleaners or abrasives on the re-usable accessories.
- Do not attempt to heat sterilize the re-usable accessories.
- Materials that are contaminated with radioactive material must be disposed of in accordance with local regulations governing radioactive materials disposal.
- Re-usable accessories shall be stored at 15-25°C and a maximum relative humidity of 85%.
- After cleaning and prior to use, inspect re-usable accessories following guidelines per Section 9.

25. IMPLANT CARD

An Implant Card (LBL011 current revision) must be given to each patient for each implantable device. Complete Implant Card per provided instructions (IFU-007 current revision).

Ensure that printed lot information, as indicated by the "LOT" symbol, matches the SIR-Spheres Y-90 resin microspheres dose lot information.

26. RADIATION SAFETY

Regulatory and local radiation usage guidelines should be followed concerning implantation and post-implantation care. A properly calibrated radiation survey meter must always be available when handling SIR-Spheres Y-90 resin microspheres to identify and control potential radioactive contamination risks.

The accessories used during the procedure are intended to provide additional protection from unintended radiation exposure for the patient and user by utilizing materials suited to use with beta radiation.

The following are sample measured thermoluminescent dosimetry (TLD) exposures to personnel.

Table 1 – Exposure Dose per Patient for Implant Preparation (Technologist)

	Trunk mSv (mrem)	Lens of the Eye mSv (mrem)	Hands mSv (mrem)
Shallow Dose (0.07mm)	0.038 (3.8)	0.12 (12)	0.32 (32)
Deep Dose (10 mm)	0.004 (0.4)	0.054 (5.4)	

Assuming handling of a 3 GBq device and dose preparation time of 30 minutes. TLDs were worn near the pelvis, on the shirt's lapel, and on the working finger.

Table 2 – Exposure Dose per Patient for Implant Procedure (Physician)

	Trunk mSv (mrem)	Lens of the Eye mSv (mrem)	Hands mSv (mrem)
Shallow Dose (0.07mm)	0.038 (3.8)	0.12 (12)	0.32 (32)
Deep Dose (10 mm)	0.004 (0.4)	0.054 (5.4)	

Assuming average patient dose of approximately 2 GBq and dose injection time of 20 minutes.

Exposure data from patients implanted with an

average of 2.1 GBq at approximately 5-6 hours post implantation at the following distances from the patient's abdomen:

Table 3 – Post-Implant Exposure

Distance	Exposure
0.25 m	18.8 µ Sv/hr
0.5 m	9.2 µ Sv/hr
1.0 m	1.5 µ Sv/hr
2.0 m	0.4 µ Sv/hr
4.0 m	<0.1 µ Sv/hr

Please note that any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established (compliance purpose)

27. CORRECTION FOR DECAY

The physical half-life of yttrium-90 is 64.1 hours. Radioactive decay factors should be applied at the time of patient dose preparation in order to calculate the true value of radioactivity present.

Table 4 – Decay Factors of SIR-Spheres Y-90 resin microspheres

Hours	Decay Factor
0.5	0.995
1	0.989
2	0.979
3	0.968
4	0.956
5	0.947
6	0.937
7	0.927
8	0.917
9	0.907
10	0.898
11	0.888
12	0.878
24	0.772

Caution: The time of the initial calibration must be converted to the user's local time.

28. INCIDENT REPORTING

Please note that any serious incident that has occurred in relation to SIR-Spheres Y-90 resin microspheres or accessories should be reported to the manufacturer (Sirtex) and the competent authority of the Member State in which the user and/or patient is established.

29. SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

The Summary of Safety and Clinical Performance (SSCP) can be consulted on the EUDAMED website:

<https://ec.europa.eu/tools/eudamed/>

30. UNIQUE DEVICE IDENTIFIER

The table below lists the UDI-DI associated to SIR Spheres® Y-90 resin microspheres device and associated Delivery System (accessories):

Table 5 – Device UDI-DI Numbers

Device name	UDI-DI
SIR-Spheres® Y-90 resin microspheres	00850014612030
Syringe Shield	00850014612085
Gen 1 Delivery Box	00850014612047
Gen 1 Delivery Set	00850014612054
Gen 1 V-Vial	00850014612061
Gen 1 V-Vial Holder	00850014612078
Gen 1 Accessories pack	00850014612092
SIROS™ Delivery Dome	00850014612009
SIROS™ Delivery Set	00850014612023
SIROS™ D-Vial Prep Set	00850014612016

SYMBOLS TABLE	
SYMBOL	SYMBOL DEFINITION
	Manufacturer
	Date of manufacture
	Consult instructions for use
	Caution
	Use-by-date
	Lot or batch code
	Catalogue number
	Serial number
	Quantity
	Ionizing radiation
	Sterilized using irradiation
	Sterilized using steam
	Sterilized using Ethylene Oxide
	Single Use Only. Indicates a medical device that is intended for use on a single patient during a single procedure.

SYMBOLS TABLE	
SYMBOL	SYMBOL DEFINITION
	Non-Sterile
	Sterile Barrier
	Protective Barrier
	Do not resterilize
	Product is not made with natural rubber latex
	Do not use if package is damaged
	Keep dry
	Temperature limit
	Relative Humidity limit
	Importer
	Authorized representative in Switzerland
	Authorized representative in the European Community
	CE mark + Notified Body identification number
	Medical Device
	Unique Device Identifier
	MR Safe
	Implant Card

SYMBOLS TABLE	
SYMBOL	SYMBOL DEFINITION
	Name of Patient
	Name and Address of Healthcare Institution / Provider
	Date of Implantation
	Informational Website for Patients
	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician or licensed healthcare practitioner

 <p>Sirtex Medical Pty Ltd Shop 6, 207 Pacific Highway St Leonards NSW 2065 Австралия Тел.: +61 2 9964 8400</p>	 <p>Sirtex Germany Manufacturing GmbH Branch Frankfurt am Main Industriepark Hoechst, building G808 65926 Frankfurt am Main Германия</p>	<p>Donawa Lifescience Consulting Srl Piazza Albania, 10 00153 Rome, Италия</p> <p>EC REP</p>	<p>MedEnvoy Switzerland Gotthardstrasse 28 6302 Zug Швейцария</p> <p>CH REP</p>	
--	---	---	---	---

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