ENGLISH

Gen 1 V-Vial Technical Description

Read the SIR-Spheres® Y-90 resin microspheres training manual before implanting SIR-Spheres Y-90 resin microspheres.

The information and steps described below are recommended for the preparation of the SIR-Spheres® Y-90 resin microspheres Instructions for Use and Training Manual before preparing SIR-Spheres Y-90 resin microspheres hereinafter referred to as SIR-Spheres.

⚠ PRECAUTIONS:

- Visually inspect all components before use, do not use if damage or wear that could impact function is observed.
- Time, distance and shielding considerations should be used to minimize exposure to radiation.
- Ensure that the distance between any puncture holes in the septum of the Gen 1 V-Vial are at least 2 mm apart.

INSTRUCTIONS

Equipment Required

The following accessories are used for SIR-Spheres dose preparation:

- Gen 1 V-Vial,
- Gen 1 V-Vial Holder, and
- Syringe Shield

The Gen 1 V-Vial (SIR-V001) is intended to provide dose preparation accessories and to hold the prescribed dose of radioactive SIR-Spheres (SIR-Y001) for implantation.

The Gen 1 V-Vial Holder (SIR-H001) is intended to provide radiation shielding for the user(s) and hold the Gen 1 V-Vial (SIR-V001) inside the Gen 1 Delivery Box (SIR-B001).

The Syringe Shield (SIR-S001) is intended to provide radiation shielding for the user(s) from SIR-Spheres contained in the syringe during dose preparation. It is used with both the Gen 1 delivery system and the SIR0S delivery system.

- Unpack SIR-Spheres, leaving the glass shipping vial in its lead pot and place it on the bench top.
- Remove the Gen 1 V-Vial from its packaging, then remove the center of the aluminium seal with forceps and wipe the rubber septum with an alcohol swab.
- Place the Gen 1 V-Vial in the acrylic Gen 1 V-Vial Holder (provided) and screw on the Gen 1 V-Vial Holder lid for stability and shielding.
- Insert a short 25 gauge needle through the rubber septum of the Gen 1 V-Vial until it just pierces the septum to create a vent. It is recommended that either a purpose-designed venting needle, or a separate filter attached to the short 25 gauge needle is used in order to prevent fluid leakage.

- Leaving the SIR-Spheres Y-90 resin microspheres glass shipping vial in the lead pot, invert the lead pot and shake vigorously before opening to re-suspend the SIR-Spheres, which will have settled during shipping.

- Quickly open the lead pot and remove the shipping vial using forceps.

- Determine the total activity of SIR-Spheres in the shipping vial using an appropriate ion chamber (dose calibrator) and then return the shipping vial to the lead pot.
- Determine the volume of SIR-Spheres that needs to be withdrawn from the shipping vial to provide the intended patient-specific activity of SIR-Spheres.
- Partially remove the aluminium seal of the SIR-Spheres shipping vial and wipe with an alcohol swab.

- Insert a 25 gauge needle through the septum of the shipping vial to create a vent, ensuring that the tip of the needle is well clear of the contents in the shipping vial. Again, it is recommended that either a purpose-designed venting needle, or a separate filter attached to the short 25 gauge needle is used in order to prevent fluid leakage.

- Attach a lubricated 21 gauge needle at least 50 mm in length to a 5 mL Luer lock syringe and place it in the acrylic Syringe Shield (provided).
- Unscrew the top of the acrylic Syringe Shield and place the syringe inside.
- Place the top of the Syringe Shield back on and tighten securely.
- Using the shielded syringe and lubricated 21 gauge needle, puncture the septum of the SIR-Spheres shipping vial and quickly draw back and forth at least six (6) times in order to re-suspend the SIR-Spheres thoroughly.
- Quickly withdraw the volume of SIR-Spheres that will provide the intended patient-specific activity, being careful of needlestick injuries, re-cap the needle safely and set aside syringe assembly in the prepared radiation work area.

- Gently agitate the shipping vial to re-disperse the microspheres and measure the activity remaining in the shipping vial with the dose calibrator; subtract the activity remaining in the shipping vial from the starting total activity in the shipping vial to determine the amount of activity that has been drawn up into the 5 mL syringe.
- If the amount of activity that has been drawn up into the 5 mL syringe is not correct, then transfer the SIR-Spheres back into the shipping vial and redraw the necessary volume of SIR-Spheres (either less or more).
- Once the correct activity has been drawn up, transfer the SIR-Spheres from the 5 mL syringe into the vented Gen 1 V-Vial in the acrylic Gen 1 V-Vial Holder. If the total volume in the 5 mL syringe is less than 3 mL, draw up enough sterile water for injection to make up to a total volume 3–5 mL before transferring the SIR-Spheres into the Gen 1 V-Vial. Ensure that the distance between any puncture holes in the septum of the Gen 1 V-Vial are at least 2 mm apart. **This step should be done ONCE only.**


- Remove the vent needle from the Gen 1 V-Vial; ensure screw on lid of the Gen 1 V-Vial Holder is secure and sit (do not force) the plug into place.
- Remove the vent needle from the shipping vial and replace the lid of the lead pot.
- The patient-specific activity of SIR-Spheres is now ready for transport to the angiography suite in which the implantation will be performed.

REUSABLE SYRINGE SHIELD CLEANING INSTRUCTIONS

- If Syringe Shield is visibly soiled, clean surfaces by spraying with bleach solution (0.5% sodium hypochlorite, NaOCl) until thoroughly wet. Allow to stand for 1 minute. Wipe with a soft cloth until all visible soil is removed.
- To disinfect, spray all surfaces with 0.5% NaOCl solution at a distance of 6-8 inches until thoroughly wet. Allow to stand for at least 1 minute. Then wipe with a soft cloth to remove any residue.
- Do not use alcohol based cleaners or abrasives on the Syringe Shield. Do not attempt to sterilize the Syringe Shield.

SYMBOLS GLOSSARY DEFINITIONS

Symbol	Definition	Symbol	Definition
	Manufacturer		Medical Device
	Date of Manufacture		Do not re-sterilize
	Consult instructions for use		Single Use Only. Indicates a medical device that is intended for use on a single patient during a single procedure.
	Caution		Product is not made with natural rubber latex
	Quantity		Do not use if package is damaged
	Catalog number		Keep dry
	Lot or batch code		Temperature limit
	Use by date		Relative Humidity limit
	Importer		Ionizing radiation
	Sterile Barrier		Unique Device Identifier
	Protective Barrier		Authorized representative in Switzerland
	Non-Sterile		Authorized representative in the European Community
	Sterilized using irradiation		CE mark + Notified Body identification number
	Sterilized using steam		MR Safe
	Serial number		Caution: Federal law (USA) restricts this device to sale by or on the order of a physician or licensed healthcare practitioner.

Manufacturer
Sirtex Medical Pty Ltd
Shop 6, 207 Pacific Highway
St Leonards, NSW 2065
Australia
Tel +612 9964 8400
Fax +612 9964 8410

Europe, Middle East & Africa
Sirtex Medical Europe GmbH
Joseph-Schumpeter-Allee 33
53227 Bonn
Germany
Tel +49 228 1840 730
Fax +49 228 1840 735

Americas
Sirtex Medical, Inc.
300 Unicorn Park Drive
Woburn, MA 01801
USA
Tel +1 888 4 SIRTEX
Tel +1 888 4 747839
Fax +1 781 721 3880

Asia Pacific
Sirtex Medical Singapore Pte Ltd
50 Science Park Road
#01-01 The Kendall
Singapore Science Park II
Singapore 117406
Tel +65 6800 6500
Fax +65 6778 9424

Donawa Lifescience
Consulting Srl
Piazza Albania, 10
00153 Rome, Italy

 
MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

 Sirtex Germany Manufacturing GmbH
Branch Frankfurt am Main
Industriepark Hoechst, building G808
65926 Frankfurt am Main
Germany

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