SIR-Spheres® Microspheres
 (Yttrium-90 Microspheres)
SIR-Y001

ENGLISH

1. DESCRIPTION

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SIR-Spheres microspheres are implanted into hepatic tumors by injection into the hepatic artery using a catheter. The SIR-Spheres microspheres distribute non-uniformly in the liver, primarily due to the unique physiological characteristics of the hepatic arterial flow, the tumor to normal liver ratio of the tissue vascularity, and the size of the tumor. The tumor usually gets higher density per unit distribution of SIR-Spheres microspheres than the normal liver. The density of SIR-Spheres microspheres in the tumor can be as high as 5 to 6 times of the normal liver tissue. In general, 1 GBq of yttrium 90/kg of tissue provides 49.67 Gy of radiation dose. Once SIR-Spheres microspheres are implanted into the liver, they are not metabolized or excreted, and they stay permanently in the liver. Each device is for single patient use.

2. INTENDED USE

SIR-Spheres Y-90 resin microspheres are intended for implantation into hepatic tumors via the hepatic artery.

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4. HOW SUPPLIED

SIR-Spheres microspheres are provided in a vial with water for injection. Each vial contains 3 GBq of $Y^{90} \pm 10\%$ (at the time of calibration) in a total of 5 cc water for injection. Each vial contains 40 – 80 million microspheres with a diameter between 20 and 60 micrometers. The vial is shipped within a 6.4 mm minimum thickness lead pot. The package consists of a crimp sealed SIR-Spheres microspheres glass vial within a lead pot, and a package insert within Type A package.

Patient specific activities are prepared in accordance with Section 11.

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These accessories are designed to shield the user during dose preparation as well as the user and patient during the implant procedure.

Sections 12 -15 provide additional safety precautions for users to ensure exposure to radiation is kept as low as reasonably achievable (ALARA).

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- had previous external beam radiation therapy to the liver
- ascites or are in clinical liver failure
- markedly abnormal liver function tests (LFTs)
- total bilirubin > 2.0 mg/dL and/or Albumin < 3.0 g/dL
- ≥ 30 Gy absorbed radiation dose to the lungs, as estimated by the Technetium-99m macroaggregated-albumin (^{99m}Tc MAA) study described in Section 10
- pre-assessment angiogram that demonstrates potential pathways for deposition of microspheres to non-target organs such as the stomach, pancreas, or bowel
- are pregnant

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Inadvertent delivery of SIR-Spheres microspheres to extra-hepatic structures such as the stomach, duodenum, gallbladder, or pancreas may result in radiation injury to these structures, including but not limited to acute abdominal pain, acute gastritis, acute cholecystitis, acute pancreatitis, and peptic ulceration. Angiographic technique must be employed to prevent the non-target delivery of SIR-Spheres microspheres to any extra-hepatic structures.

• Radioembolization Induced Liver Disease (REILD)

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- Reduced liver functional reserve due to steatosis, steatohepatitis, hepatitis or cirrhosis
- Elevated baseline bilirubin level
- Prior hepatic resection
- Prior liver directed therapy
- Extensive prior treatment with systemic chemotherapy and/or biologic therapies

• Radiation Pneumonitis

High levels of implanted radiation and/or excessive shunting to the lung may lead to radiation pneumonitis. The lung radiation dose must be limited to ≤ 30 Gy for a single treatment session and ≤ 50 Gy cumulative dose.

• Other Tumor Types

The safety, performance, and benefit/risk profile of SIR-Spheres in treating particular types of tumors outside the indications for use has yet to be established.

8. PRECAUTIONS

- Safety and effectiveness of this device in pregnant women, nursing mothers or children have not been established.
- A SPECT or PET scan of the upper abdomen be performed immediately after

implantation of SIR-Spheres microspheres. The SPECT or PET scan will detect radiation from the yttrium-90 to confirm placement of the microspheres in the liver.

- This product is radioactive. The use of this device is regulated, and local regulations must be followed when handling this device.
- Accepted radiation protection techniques should be used to protect staff when handling both the isotope and the patient.
- Patients may experience gastric problems following treatment, but proton pump inhibitors (PPI) or histamine H2-receptor antagonists (H-2 blocking agents) may be used the day before implantation of SIR-Spheres microspheres and continued as needed to reduce gastric complications.
- Patients may experience abdominal pain immediately after administration of SIR-Spheres microspheres and pain relief may be required.
- SIR-Spheres microspheres demonstrated a mild sensitization potential when tested dermally in an animal model.

9. ADVERSE EVENTS

When the patient is treated with proper technique, without excessive radiation to any organ, the common adverse events after receiving the SIR-Spheres microspheres are fever, transient decrease of hemoglobin, transient thrombocytopenia, mild to moderate abnormality of liver function tests (mild increase in aspartate aminotransferase, alkaline phosphatase, bilirubin), abdominal pain, nausea, vomiting, and diarrhea.

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Radiation induced liver disease

(REILD): REILD is a rare complication following Selective Internal Radiation Therapy (SIRT). REILD is characterized by a well-defined constellation of temporal, clinical, biochemical, and histopathologic findings. It typically manifests approximately 4 to 8 weeks post-SIRT and is characterized clinically by jaundice and ascites in the absence of tumor progression or bile duct obstruction.

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In the treatment of REILD, low molecular weight heparin may also be considered but both corticosteroids and heparin may only be useful if commenced very early in the course of the disease. See also Section 7 Warnings.

¹ Gil-Alzugaray et al. Hepatology, Vol 57, No. 3, 2013.

10. PATIENT SELECTION AND PRETREATMENT TESTING

The potential benefit of selective internal radiation therapy for disease control within the liver is realized in patients with appropriate vascular anatomy of the liver and surrounding tissue, relatively good liver function or functional liver reserve, low lung shunt, and appropriate activity prescription.

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- An hepatic angiogram to establish arterial anatomy of the liver
- ^{99m}Tc MAA Scan to determine the percent lung shunting to evaluate lung radiation dose exposure
- Biochemical tests of liver function
- CT, PET/CT, or MRI to determine the extent of the disease, lung, tumor and liver volumes for dosimetry calculations

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10.3 Lung Shunt Calculation Procedure

Calculate the lung shunt fraction (L) using following formula:

$$L = \left(\frac{\text{Lung Counts}}{\text{Liver Counts} + \text{Lung Counts}} \right)$$

Equation 1

In order to optimize the risk versus benefit for patients receiving SIR-Spheres microspheres, limiting the radiation exposure to the lungs to ≤ 30 Gy is required. The calculation of estimated radiation exposure to the lungs is given by the following formula:

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The approach of lobar treatment versus whole liver treatment with SIR-Spheres microspheres is based on the presence of visible tumors on pre-treatment CT or MR imaging. If liver tumors are only visible in one lobe, then SIR-Spheres microspheres should be administered to that lobe only, thus sparing the contralateral lobe from unnecessary internal radiation.

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Equation 4

Where:

BSA = Body Surface Area
 H = height in meters
 W = weight in kilograms

11.1.1 BSA prescribed activity calculation for whole liver / bilobar treatment

$$A_{Admin} = (BSA - 0.2) + \left(\frac{V_{Tumor}}{V_{Tumor} + V_{Normal\ liver}} \right)$$

Equation 5

Where:

A_{Admin} [GBq] = SIR-Spheres activity to implant
 V_{Tumor} = volume of tumor
 $V_{Normal\ Liver}$ = volume of non-tumor tissue in the treated volume

11.1.2 BSA prescribed activity calculation for lobar or super-selective treatment

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Equation 6

Where:

A_{Admin} = SIR-Spheres activity to implant [GBq]
 V_{Tumor} = volume of tumor in the portion of the liver being treated (i.e., lobe)
 $V_{Treated}$ = volume of the portion of the liver being treated AND the tumor volume
 V_{Liver} = total volume of the whole liver, including tumor
 BSA = Body Surface Area as per Equation 4

11.2 Partition Model

This method involves selecting safe radiation doses to the normal liver and lung and implanting the maximum activity that will not exceed these limits. The dose received by the tumor has no upper limit.

The partition model should be used where the tumor mass is a discrete area within the liver. The technique requires two measurements to be made:

1. Measurement of the volume of tumor and normal perfused liver determined from a CT or MR scan
2. Measurement of the proportions of ^{99m}Tc MAA activity that lodges in the tumor, normal perfused liver and lung as determined from imaging.

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- Determine the volumes of the normal perfused liver and tumor and convert each volume to mass²
- Determine the volume of the lung and convert to mass³
- Using the ^{99m}Tc MAA scan, determine the activity in the lung, tumor and perfused normal liver

- Determine the Tumor to Normal Activity Ratio (TNR) calculated as activity per unit mass of organ or tissue using the following equation:

$$TNR = \frac{A_{Tumor}/M_{Tumor}}{A_{PNL}/M_{PNL}}$$

$$= \frac{\text{average counts or activity in tumor}}{\text{average counts or activity in perfused normal liver}}$$

Equation 7

Where:

A_{Tumor} = Activity in tumor
 M_{Tumor} = Mass of tumor
 A_{PNL} = Activity in the perfused normal liver
 M_{PNL} = Mass of the perfused normal liver

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Equation 8

Where:

A_{Admin} = SIR-Spheres activity to implant [GBq]
 D_{PNL} = Mean absorbed dose to perfused normal liver [Gy]
 M_{Tumor} = Mass of tumor [g]
 M_{PNL} = Mass of perfused normal liver [g]
 L = Lung Shunt Fraction

and,

$$D_{Tumor} = TNR * D_{PNL}$$

Equation 9

Where:

D_{Tumor} = Desired absorbed dose to tumor [Gy]

12. DOSE PREPARATION PROCEDURE

12.1 Accessories provided by Sirtex:

- V-Vial and capped 21 G needle (SIR-V001)
- V-Vial Holder (SIR-H001)
- Syringe Shield (SIR-S001)

12.2 Additional supplies needed:

- Sterile water for Injection
- Two 25 G vial needles with filters
- 5 mL syringe
- Alcohol swabs
- Forceps
- Dose Calibrator (ion chamber)
- Radiation survey meter
- Radiation spill kit

12.3 Dose Prep Procedure

- a. unpack SIR-Spheres microspheres, leaving the glass shipping vial in its lead pot and place on the benchtop.
- b. Remove the center of the aluminium seal from the sterile V-Vial with forceps and wipe the septum with an alcohol swab.
- c. Place the V-Vial in the V-Vial holder.
- d. Insert a short 25-gauge needle with filter through the septum of the V-Vial until it just pierces the septum to create a vent.
- e. Leaving the SIR-Spheres microspheres shipping vial in the lead pot and shake to re-suspend the SIR-Spheres microspheres. Resuspension ensures a homogenous solution for dose preparation.
- f. Open the lead pot and remove the shipping vial using forceps.
- g. Determine the total activity of SIR-Spheres microspheres in the shipping vial using a dose calibrator, then return the shipping vial to the lead pot.

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- h. Determine the volume to be withdrawn to provide the intended patient specific activity.
- i. Remove the center of the aluminium seal of the SIR-Spheres microspheres shipping vial, wipe septum with alcohol swab.
- j. Insert a 25-gauge needle with filter through the septum of the shipping vial to create a vent, ensuring the needle is well clear of the contents in the shipping vial.
- k. Use a shielded syringe with a 21-gauge needle at least 70 mm long to puncture the septum of the SIR-Spheres microspheres shipping vial, and quickly draw back and forth several times in order to mix the SIR-Spheres microspheres thoroughly.
- l. Withdraw the precalculated volume.
- m. Verify the patient dose in the V-Vial by remeasuring the activity in the shipping vial with the dose calibrator, and correct, if necessary.
- n. Transfer the patient specific activity into the vented V-Vial in V-Vial holder.
- o. Remove the vent and place the black plug securely into the opening.

The patient specific activity is now ready for transport to the SIR-Spheres microspheres implantation room.

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13.2 Additional supplies needed:

- Two 20 mL Luer Lock Syringes filled with non-ionic solution (either 5% glucose/ dextrose or water for injection).
 - Do not use saline
 - Patients with diabetes should be monitored for hyperglycemia if 5% dextrose/glucose is used to implant SIR-Spheres
- One 20 mL Luer Lock Syringe filled with non-ionic contrast
- Microcatheter
 - The microcatheters should have an inner diameter of at least 0.021" and a 45° tip configuration for vessels with a 90° angle of origin.
- Sterile absorbent material for protection of angiography suite
- Alcohol swabs
- Forceps or hemostat
- Radiation survey meter
- Radiation spill kit

13.3 Implant Procedure

The hepatic artery catheter is inserted by a trained interventional radiologist under x-ray guidance. This method allows complete control of exactly where the catheter is placed and allows routine checking of the catheter position throughout the implant procedure.

It is essential that SIR-Spheres microspheres are not delivered to other organs, in particular the pancreas, stomach, or duodenum. If there is any possibility of SIR-Spheres microspheres passing down the Gastroduodenal artery (GDA) then the implantation must not proceed. It may be preferable to block shunting vessels with an intraluminal coil or other agent to prevent SIR-Spheres microspheres from flowing to non-target organs.

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SIR-Spheres microspheres must be delivered slowly at a rate of no more than 5 ml per minute. Rapid delivery may result in a more concentrated suspension of SIR-Spheres that can cause clogging in the microcatheter or reflux back down the hepatic artery and into other organs. At the conclusion of the procedure, the catheter is removed.

14. DISASSEMBLY

- Once the infusion is complete, remove lid from the delivery box.
- Do not disconnect the patient catheter from the tubing set.
- Carefully remove catheter from the patient, while still attached to the Delivery Set. Use care when removing and handling the catheter as it may be contaminated with radioactivity. Coil the catheter and wrap within a sterile towel.
- Leaving everything attached, use a hemostat to push down the needles into the shielded V-Vial and place it (along with all potentially contaminated procedural accessories) into the designated receptacle for radioactive material.
- Radioactive materials must be stored in accordance with local regulations governing radioactive materials storage.

15. CLEANUP AND WASTE DISPOSAL

- Following dose preparation and delivery, survey equipment and people using appropriate radiation survey meter to ensure any contamination is contained.
- If re-usable accessories (i.e., Syringe Shield or Delivery box) are visibly soiled, clean surfaces by spraying with cleaning solution (0.5% sodium hypochlorite, NaOCl) until thoroughly wet. Allow it to stand for 1 minute. Wipe with a soft cloth until all visible soil is removed.
- To disinfect, spray all surfaces with 0.5% NaOCl solution at 6-8 inches until thoroughly wet. Allow to stand for at least 1 minute. Then wipe with a soft cloth to remove any residue.
- Do not use alcohol-based cleaners or abrasives on the reusable accessories.
- Do not attempt to heat sterilize the reusable accessories.
- Materials that are contaminated with radioactive material must be disposed of in accordance with local regulations governing radioactive materials disposal.

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Table 2 – Exposure Dose per Patient for Implant Procedure (Physician)

	Trunk mSv (mrem)	Lens of the Eye mSv (mrem)	Hands mSv (mrem)
Shallow Dose (0.07 mm)	0.038 (3.8)	0.12 (12)	0.32 (32)
Deep Dose (10 mm)	0.004 (0.4)	0.054 (5.4)	

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SYMBOLS TABLE	
SYMBOL	SYMBOL DEFINITION
	Manufacturer
	Date of manufacture
	Consult instructions for use
	Caution
	Use-by-date
	Lot or batch code
	Catalog number
	Serial number
	Quantity
	Sterilized using irradiation
	Sterilized using steam

SYMBOL	SYMBOL DEFINITION
	Ionizing radiation
	Single Use Only. Indicates a medical device that is intended for use on a single patient during a single procedure
	Do not resterilize
	Product is not made with natural rubber latex
	Do not use if package is damaged
	Keep dry
	Temperature limit
	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician or licensed healthcare practitioner
	Importer
	Authorized representative in Switzerland
	Authorized representative in the European Community
	CE mark + Notified Body identification number

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ENGLISH (AUSTRALIA)

1. DESCRIPTION

SIR-Spheres microspheres consist of biocompatible microspheres containing yttrium-90. Yttrium-90 is a high-energy pure beta-emitting isotope with no primary gamma emission. The maximum energy of the beta particles is 2.27 MeV with a mean of 0.93 MeV. The half-life is 64.1 hours. The maximum range of emissions in tissue is 11 mm with a mean of 2.5 mm.

SIR-Spheres microspheres are implanted into hepatic tumours by injection into the hepatic artery using a catheter. The SIR-Spheres microspheres distribute non-uniformly in the liver, primarily due to the unique physiological characteristics of the hepatic arterial flow, the tumour-to-normal-liver ratio of the tissue vascularity, and the size of the tumour. The tumour usually gets higher density-per-unit distribution of SIR-Spheres microspheres than the normal liver. The density of SIR-Spheres microspheres in the tumour can be as high as 5 to 6 times that of the normal liver tissue. In general, 1 GBq of yttrium 90 / kg of tissue provides 49.67 Gy of radiation dose. Once SIR-Spheres microspheres are implanted into the liver, they are not metabolised or excreted, and they stay permanently in the liver. Each device is for single patient use.

2. INTENDED USE

SIR-Spheres Y-90 resin microspheres are intended for implantation into hepatic tumours via the hepatic artery.

3. INDICATIONS FOR USE

SIR-Spheres Y-90 resin microspheres are indicated for the treatment of patients with advanced non-operable liver cancer.

4. HOW SUPPLIED

SIR-Spheres microspheres are provided in a vial with water for injection. Each vial contains 3 GBq of $Y^{90} \pm 10\%$ (at the time of calibration) in a total of 5 cc water for injection. Each vial contains 40 – 80 million microspheres with a diameter between 20 and 60 micrometres. The vial is shipped within a 6.4 mm minimum thickness lead pot. The package consists of a crimp-sealed SIR-Spheres microspheres glass vial within a lead pot, and a package insert within Type A package.

Patient-specific activities are prepared in accordance with Section 11.

The vial and its contents should be stored inside its transportation container at room temperature (15-25 °C; 59-77 °F).

The calibration date (for radioactive contents) and the expiration information are printed on the vial label. The useful life of the SIR-Spheres microspheres ends 24 hours after the time of calibration.

5. ACCESSORIES

Accessories that may be used for the dose preparation and implant procedure include:

- Delivery Set (SIR-D001)
- V-Vial and capped 21 G needle (SIR-V001)
- V-Vial Holder (SIR-H001)
- Delivery Box (SIR-B001)
- Syringe Shield (SIR-S001)

These accessories are designed to shield the user during dose preparation as well as the user and patient during the implant procedure. Sections 12-15 provide additional safety precautions for users to ensure exposure to radiation is kept as low as reasonably achievable (ALARA).

6. CONTRAINDICATIONS

SIR-Spheres microspheres are contraindicated in patients who have:

- had previous external beam radiation therapy to the liver
- ascites or are in clinical liver failure
- markedly abnormal liver function tests (LFTs)
- total bilirubin > 2.0 mg/dL and/or Albumin < 3.0 g/dL
- ≥ 30 Gy absorbed radiation dose to the lungs, as estimated by the Technetium-99m macroaggregated-albumin (^{99m}Tc MAA) study described in Section 10
- pre-assessment angiogram that demonstrates potential pathways for deposition of microspheres to non-target organs such as the stomach, pancreas or bowel
- are pregnant

7. WARNINGS

• Non-Target Delivery of SIR-Spheres microspheres

Inadvertent delivery of SIR-Spheres microspheres to extra-hepatic structures such as the stomach, duodenum, gallbladder or pancreas may result in radiation injury to these structures, including but not limited to acute abdominal pain, acute gastritis, acute cholecystitis, acute pancreatitis and peptic ulceration. Angiographic technique must be employed to prevent the non-target delivery of SIR-Spheres microspheres to any extra-hepatic structures.

• Radioembolisation-Induced Liver Disease (REILD)

Delivery of excessive radiation to the normal liver parenchyma may result in REILD. The risk of REILD may also be increased in patients with pre-existing liver disease. Consideration should be given to reducing the prescribed activity of SIR-Spheres microspheres in the following clinical settings¹:

- Reduced liver functional reserve due to steatosis, steatohepatitis, hepatitis or cirrhosis
- Elevated baseline bilirubin level
- Prior hepatic resection
- Prior liver-directed therapy
- Extensive prior treatment with systemic chemotherapy and/or biologic therapies
- **Radiation Pneumonitis**
High levels of implanted radiation and/or excessive shunting to the lung may lead to radiation pneumonitis. The lung radiation dose must be limited to ≤ 30 Gy for a single treatment session and ≤ 50 Gy cumulative dose.
- **Other Tumour Types**
The safety, performance and benefit/risk profile of SIR-Spheres in treating particular types of tumours outside the indications for use has yet to be established.

8. PRECAUTIONS

- Safety and effectiveness of this device in pregnant women, nursing mothers or children have not been established.
- A SPECT or PET scan of the upper abdomen should be performed immediately after implantation of SIR-Spheres

microspheres. The SPECT or PET scan will detect radiation from the yttrium-90 to confirm placement of the microspheres in the liver.

- This product is radioactive. The use of this device is regulated, and local regulations must be followed when handling this device.
- Accepted radiation protection techniques should be used to protect staff when handling both the isotope and the patient.
- Patients may experience gastric problems following treatment, but proton pump inhibitors (PPI) or histamine H2-receptor antagonists (H-2 blocking agents) may be used the day before implantation of SIR-Spheres microspheres and continued as needed to reduce gastric complications.
- Patients may experience abdominal pain immediately after administration of SIR-Spheres microspheres and pain relief may be required.
- SIR-Spheres microspheres demonstrated a mild sensitisation potential when tested dermally in an animal model.

9. ADVERSE EVENTS

When the patient is treated with proper technique, without excessive radiation to any organ, the common adverse events after receiving the SIR-Spheres microspheres are fever, transient decrease of haemoglobin, transient thrombocytopenia, mild to moderate abnormality of liver function tests (mild increase in aspartate aminotransferase, alkaline phosphatase, bilirubin), abdominal pain, nausea, vomiting and diarrhoea.

Potential Serious Adverse Events Due to High Radiation

Acute Pancreatitis: causes immediate severe abdominal pain. Verify by SPECT or PET of the abdomen and test for serum amylase.

Radiation Pneumonitis: causes excessive non-productive cough. Verify by X-ray evidence of pneumonitis.

Acute Gastritis: causes abdominal pain. Verify by standard methods to diagnose gastric ulceration.

Acute Cholecystitis: causes significant upper abdominal pain and may require cholecystectomy for resolution. Verify by appropriate imaging studies.

Radioembolisation-induced liver disease (REILD): REILD is a rare complication following Selective Internal Radiation Therapy (SIRT). REILD is characterised by a well-defined constellation of temporal, clinical, biochemical and histopathologic findings. It typically manifests approximately 4 to 8 weeks post-SIRT and is characterised clinically by jaundice and ascites in the absence of tumour progression or bile duct obstruction.

The typical biochemical picture of REILD is an elevated bilirubin (> 3 mg/dL) in almost all cases, elevated alkaline phosphatase (ALP) and gamma-glutamyl transpeptidase (GGT) in most cases, accompanied by virtually no change in the transaminases (AST and ALT). If liver biopsy is performed, the typical histological appearance is of sinusoidal obstruction that may resemble veno-occlusive disease. REILD may occur in both non-cirrhotic and cirrhotic patients.

Prophylactic treatment with methylprednisolone and ursodeoxycholic acid starting on the day of SIRT and continued for two months may reduce the incidence of REILD.

In the treatment of REILD, low-molecular-weight heparin may also be considered but both corticosteroids and heparin may only be useful if commenced very early in the course of the disease. See also Section 7 Warnings.

¹ Gil-Alzugaray et al. Hepatology, Vol 57, No. 3, 2013.

10. PATIENT SELECTION AND PRE-TREATMENT TESTING

The potential benefit of selective internal radiation therapy for disease control within the liver is realised in patients with appropriate vascular anatomy of the liver and surrounding tissue, relatively good liver function or functional liver reserve, low lung shunt and appropriate activity prescription.

10.1 Patient Tests Before Treatment with SIR-Spheres microspheres

The following tests are performed before treatment to ensure the patient is suitable for treatment with SIR-Spheres Y-90 resin microspheres and determine the appropriate dose of SIR-Spheres to prescribe:

- A hepatic angiogram to establish arterial anatomy of the liver
- ^{99m}Tc MAA Scan to determine the percentage of lung shunting to evaluate lung radiation dose exposure
- Biochemical tests of liver function
- CT, PET/CT or MRI to determine the extent of the disease, lung, tumour and liver volumes for dosimetry calculations

10.2 Technique for Performing the Intra-Hepatic Technetium MAA Scan

To assess arterial perfusion of the liver and the fraction of radiopharmaceutical tracer that will pass through the liver and lodge in the lungs:

Inject about 150 MBq of ^{99m}Tc MAA into the hepatic artery via a catheter.

Use a large FOV gamma camera and obtain images of the thorax and abdomen (with the same acquisition time).

Draw Region of Interest (ROI) around the whole liver and the whole lung and get the total counts for the lung and the liver.

10.3 Lung Shunt Calculation Procedure

Calculate the lung shunt fraction (*L*) using following formula:

$$L = \left(\frac{\text{Lung Counts}}{\text{Liver Counts} + \text{Lung Counts}} \right)$$

Equation 1

In order to optimise the risk versus benefit for patients receiving SIR-Spheres microspheres, limiting the radiation exposure to the lungs to ≤ 30 Gy is required. The calculation of estimated radiation exposure to the lungs is given by the following formula:

Activity that may potentially reach the lung:

$$A_{\text{lung}} = A_{\text{total}} \times L$$

Equation 2

Where:

A_{lung} = lung activity [GBq]
 A_{total} = total prescribed activity [GBq]
 L = lung shunt fraction

The resulting lung dose, given that a given amount of activity shunts from the liver to the lung:

$$D_{\text{lung}} = \frac{49670 \times A_{\text{lung}}}{M_{\text{lung}}}$$

Equation 3

Where:

D_{lung} = lung dose [Gy]
 A_{lung} = lung activity [GBq]
 M_{lung} = mass of the lung [g]

11. CALCULATION OF INDIVIDUAL DOSE

There are two accepted methods for calculating the patient radiation dose, these being the Body Surface Area (BSA) model and the partition model.

11.1 The Body Surface Area (BSA) Method

The BSA method varies yttrium-90 activity according to the size of the patient and the size of the tumour within the liver. The BSA method can be used when treating the volume of a single treated lobe, as well as that of the entire liver.

The approach of lobar treatment versus whole-liver treatment with SIR-Spheres microspheres is based on the presence of visible tumours on pre-treatment CT or MR imaging. If liver tumours are only visible in one lobe, then SIR-Spheres microspheres should be administered to that lobe only, thus sparing the contralateral lobe from unnecessary internal radiation.

BSA must first be determined and is calculated from the following equation:

$$BSA = 0.20247 \times H^{0.725} \times W^{0.425}$$

Equation 4

Where:

BSA = Body Surface Area
 H = height in metres
 W = weight in kilograms

11.1.1 BSA prescribed-activity calculation for whole-liver / bilobar treatment

$$A_{\text{Admin}} = (BSA - 0.2) + \left(\frac{V_{\text{Tumour}}}{V_{\text{Tumour}} + V_{\text{Normal Liver}}} \right)$$

Equation 5

Where:

A_{Admin} [GBq] = SIR-Spheres activity to implant
 V_{Tumour} = volume of tumour
 $V_{\text{Normal Liver}}$ = volume of non-tumour tissue in the treated volume

11.1.2 BSA prescribed-activity calculation for lobar or super-selective treatment

In patients who receive lobar or segmental treatment with SIR-Spheres microspheres, the prescribed activity must be reduced in accordance with the size of the portion of the liver being treated.

$$A_{\text{Admin}} = \left[BSA - 0.2 + \left(\frac{V_{\text{Tumour}}}{V_{\text{Treated}}} \right) \right] \times \left[\frac{V_{\text{Treated}}}{V_{\text{Liver}}} \right]$$

Equation 6

Where:

A_{Admin} = SIR-Spheres activity to implant [GBq]
 V_{Tumour} = volume of tumour in treated volume (i.e. lobe)
 V_{Treated} = volume of non-tumour tissue in the treated volume (i.e. lobe)

V_{Liver} = total volume of the whole liver, including tumour
 BSA = Body Surface Area as per Equation 5

11.2 Partition Model

This method involves selecting safe radiation doses to the normal liver and lung and implanting the maximum activity that will not exceed these limits. The dose received by the tumour has no upper limit.

The partition model should be used where the tumour mass is a discrete area within the liver. The technique requires two measurements to be made:

1. Measurement of the volume of tumour and normal perfused liver determined from a CT or MR scan
2. Measurement of the proportions of ^{99m}Tc MAA activity that lodges in the tumour, normal perfused liver and lung as determined from imaging.

Therefore, to calculate the activity to be implanted, it is necessary to:

- Determine the volumes of the normal perfused liver and tumour and convert each volume to mass²
- Determine the volume of the lung and convert to mass³

- Using the ^{99m}Tc MAA scan, determine the activity in the lung, tumour and perfused normal liver
- Determine the Tumour-to-Normal-Activity Ratio (TNR) calculated as activity per unit mass of organ or tissue using the following equation:

$$TNR = \frac{A_{\text{Tumour}}/M_{\text{Tumour}}}{A_{\text{PNL}}/M_{\text{PNL}}}$$

$$= \frac{\text{average counts or activity in tumour}}{\text{average counts or activity in perfused normal liver}}$$

Equation 7

Where:

A_{Tumour} = Activity in tumour
 M_{Tumour} = Mass of tumour
 A_{PNL} = Activity in the perfused normal liver
 M_{PNL} = Mass of the perfused normal liver

To calculate the total activity to be implanted, use the following equations. The activity required should be calculated to accommodate the lung and normal tissue doses as the limiting factors.

$$A_{\text{Admin}} = \frac{D_{\text{PNL}} (TNR * M_{\text{Tumour}} + M_{\text{PNL}})}{49670 * (1 - L)}$$

Equation 8

Where:

A_{Admin} = SIR-Spheres activity to implant [GBq]
 D_{PNL} = Mean absorbed dose to perfused normal liver [Gy]
 M_{Tumour} = Mass of tumour [g]
 M_{PNL} = Mass of perfused normal liver [g]
 L = Lung Shunt Fraction

and,

$$D_{\text{Tumour}} = TNR * D_{\text{PNL}}$$

Equation 9

Where:

D_{Tumour} = Desired absorbed dose to tumour [Gy]

12. DOSE PREPARATION PROCEDURE

12.1 Accessories provided by Sirtex:

- V-Vial and capped 21 G needle (SIR-V001)
- V-Vial Holder (SIR-H001)
- Syringe Shield (SIR-S001)

12.2 Additional supplies needed:

- Sterile water for Injection
- Two 25 G vent needles with filters
- 5 mL syringe
- Alcohol swabs
- Forceps
- Dose Calibrator (ion chamber)
- Radiation survey meter
- Radiation spill kit

12.3 Dose Prep Procedure

- a. unpack SIR-Spheres microspheres, leaving the glass shipping vial in its lead pot and place on the benchtop.
- b. Remove the centre of the aluminium seal from the sterile v-vial with forceps and wipe the septum with an alcohol swab.
- c. Place the V-Vial in the V-Vial holder.
- d. Insert a short 25-gauge needle with filter through the septum of the v-vial until it just pierces the septum to create a vent.
- e. Leaving the SIR-Spheres microspheres shipping vial in the lead pot and shake to re-suspend the SIR-Spheres microspheres. Resuspension ensures a homogeneous solution for dose preparation.
- f. Open the lead pot and remove the shipping vial using forceps.

² ICRP, 2015. Radiation Dose to Patients from Radiopharmaceuticals: A Compendium of Current Information Related to Frequently Used Substances. ICRP Publication 128. Ann. ICRP 44(2S).

³ Kappadath SC, Lopez BP, Salem R, Lam MG. Lung shunt and lung dose calculation methods for radioembolization treatment planning. Q J Nucl Med Mol Imaging. 2021 Mar;65(1):32-42.

- g. Determine the total activity of SIR-Spheres microspheres in the shipping vial using a dose calibrator, then return the shipping vial to the lead pot.
- h. Determine the volume to be withdrawn to provide the intended patient-specific activity.
- i. Remove the centre of the aluminium seal of the SIR-Spheres microspheres shipping vial, wipe septum with alcohol swab.
- j. Insert a 25-gauge needle with filter through the septum of the shipping vial to create a vent, ensuring the needle is well clear of the contents in the shipping vial.
- k. Use a shielded syringe with a 21-gauge needle at least 70 mm long to puncture the septum of the SIR-Spheres microspheres shipping vial, and quickly draw back and forth several times in order to mix the SIR-Spheres microspheres thoroughly.
- l. Withdraw the pre-calculated volume.
- m. Verify the patient dose in the V-Vial by remeasuring the activity in the shipping vial with the dose calibrator, and correct, if necessary.
- n. Transfer the patient-specific activity into the vented V-Vial in V-Vial holder.
- o. Remove the vent and place the black plug securely into the opening.

The patient-specific activity is now ready for transport to the SIR-Spheres microspheres implantation room.

13. IMPLANT PROCEDURE

Doctors must refer to the Sirtex Medical Pty Ltd Training Manual for delivering SIR-Spheres resin microspheres before attempting to implant this device.

13.1 Accessories provided by Sirtex:

- Delivery Set (SIR-D001)
- Delivery Box (SIR-B001)
- V-Vial (SIR-V001)
- V-Vial Holder (SIR-H001)

13.2 Additional supplies needed:

- Two 20 mL Luer Lock Syringes filled with non-ionic solution (either 5% glucose/dextrose or water for injection).
 - Do not use saline
 - Patients with diabetes should be monitored for hyperglycaemia if 5% dextrose/glucose is used to implant SIR-Spheres
- One 20 mL Luer Lock Syringe filled with non-ionic contrast
- Microcatheter
 - The microcatheters should have an inner diameter of at least 0.021" and a 45° tip configuration for vessels with a 90° angle of origin.
- Sterile absorbent material for protection of angiography suite
- Alcohol swabs
- Forceps or hemostat
- Radiation survey meter
- Radiation spill kit

13.3 Implant Procedure

The hepatic artery catheter is inserted by a trained interventional radiologist under X-ray guidance. This method allows complete control of exactly where the catheter is placed and allows routine checking of the catheter position throughout the implant procedure.

It is essential that SIR-Spheres microspheres are not delivered to other organs, in particular the pancreas, stomach or duodenum. If there is any possibility of SIR-Spheres microspheres passing down the gastroduodenal artery (GDA), then the implantation must not proceed. It may be preferable to block shunting vessels with an intraluminal coil or other agent to prevent SIR-Spheres microspheres from flowing to non-target organs.

Note: Virtually all complications from SIR-Spheres microspheres arise from the inadvertent delivery of SIR-Spheres microspheres into small blood vessels that go to the pancreas, stomach or duodenum.

The radiologist must repeatedly check the position of the catheter during the procedure to ensure it remains correctly sited and that reflux of the SIR-Spheres microspheres into other organs does not occur. This is performed by injecting contrast medium through the B line port of the SIR-Spheres Delivery Set during the delivery of SIR-Spheres microspheres. Contrast must not be administered on the D line port.

SIR-Spheres microspheres must be delivered slowly at a rate of no more than 5 mL per minute. Rapid delivery may result in a more concentrated suspension of SIR-Spheres, which can cause clogging in the microcatheter or reflux back down the hepatic artery and into other organs. At the conclusion of the procedure, the catheter is removed.

14. DISASSEMBLY

- Once the infusion is complete, remove lid from the delivery box.
- Do not disconnect the patient catheter from the tubing set.
- Carefully remove the catheter from the patient, while still attached to the Delivery Set. Use care when removing and handling the catheter as it may be contaminated with radioactivity. Coil the catheter and wrap within a sterile towel.
- Leaving everything attached, use a hemostat to push down the needles into the shielded V-Vial and place it (along with all potentially contaminated procedural accessories) into the designated receptacle for radioactive material.
- Radioactive materials must be stored in accordance with local regulations governing radioactive materials storage.

15. CLEAN-UP AND WASTE DISPOSAL

- Following dose preparation and delivery, survey equipment and people using appropriate radiation survey meter to ensure any contamination is contained.
- If re-usable accessories (i.e. Syringe Shield or Delivery box) are visibly soiled, clean surfaces by spraying with cleaning solution (0.5% sodium hypochlorite, NaOCl) until thoroughly wet. Allow it to stand for 1 minute. Wipe with a soft cloth until all visible soil is removed.
- To disinfect, spray all surfaces with 0.5% NaOCl solution at 6-8 inches until thoroughly wet. Allow to stand for at least 1 minute. Then wipe with a soft cloth to remove any residue.
- Do not use alcohol-based cleaners or abrasives on the reusable accessories.
- Do not attempt to heatsterilise the reusable accessories.
- Materials that are contaminated with radioactive material must be disposed of in accordance with local regulations governing radioactive materials disposal.

16. RADIATION SAFETY

Regulatory and local radiation usage guidelines should be followed concerning implantation and post-implantation care. A properly calibrated radiation survey meter must always be available when handling SIR-Spheres to identify and control potential radioactive contamination risks.

The following are sample measured thermoluminescent dosimetry (TLD) exposures to personnel.

Table 1 – Exposure Dose per Patient for Implant Preparation (Technologist)

	Trunk mSv (mrem)	Lens of the Eye mSv (mrem)	Hands mSv (mrem)
Shallow Dose (0.07 mm)	0.027 (2.7)	0.026 (2.6)	0.35 (35)
Deep Dose (10 mm)	0.003 (0.3)	0.004 (0.4)	

Assuming handling of a 3 GBq device and dose preparation time of 30 minutes. TLDs were worn near the pelvis, on the shirt's lapel, and on the working finger.

Table 2 – Exposure Dose per Patient for Implant Procedure (Doctor)

	Trunk mSv (mrem)	Lens of the Eye mSv (mrem)	Hands mSv (mrem)
Shallow Dose (0.07 mm)	0.038 (3.8)	0.12 (12)	0.32 (32)
Deep Dose (10 mm)	0.004 (0.4)	0.054 (5.4)	

Assuming average patient dose of approximately 2 GBq and dose injection time of 20 minutes.

Exposure data from patients implanted with an average of 2.1 GBq at approximately 5-6 hours post implantation at the following distances from the patient's abdomen:

Table 3 – Post-Implant Exposure

Distance	Exposure
0.25 m	18.8 µSv/hr
0.5 m	9.2 µSv/hr
1.0 m	1.5 µSv/hr
2.0 m	0.4 µSv/hr
4.0 m	< 0.1 µSv/hr

17. CORRECTION FOR DECAY

The physical half-life of yttrium-90 is 64.1 hours. Radioactive decay factors should be applied at the time of patient dose preparation in order to calculate the true value of radioactivity present.


Table 4 – Decay Factors of SIR-Spheres microspheres

Hours	Decay Factor
0.5	0.995
1	0.989
2	0.979
3	0.968
4	0.956
5	0.947
6	0.937
7	0.927
8	0.917
9	0.907
10	0.898
11	0.888
12	0.878
24	0.772



Caution: The time of the initial calibration must be converted to the user's local time.

CH REP

MedEnvoy Switzerland
 Gotthardstrasse 28
 6302 Zug
 Switzerland

SYMBOLS TABLE	
SYMBOL	SYMBOL DEFINITION
	Manufacturer
	Date of manufacture
	Consult instructions for use
	Caution
	Use-by date
	Lot or batch code
	Catalogue number
	Serial number
	Quantity
	Sterilised using irradiation
	Sterilised using steam

SYMBOL	SYMBOL DEFINITION
	Ionising radiation
	Single Use Only. Indicates a medical device that is intended for use on a single patient during a single procedure
	Do not resterilise
	Product is not made with natural rubber latex
	Do not use if package is damaged
	Keep dry
	Temperature limit
	Caution: Australian law restricts this device to sale by or on the order of a doctor or licensed healthcare practitioner
	Importer
	Authorised representative in Switzerland
	Authorised representative in the European Community
	CE mark + Notified Body identification number

IMPLANT CARD SYMBOLS	
SYMBOL	SYMBOL DEFINITION
	Implant Card
	Patient name or patient ID
	Date of implantation
	Name and Address of the implanting healthcare institution/provider
	Information website for patients
	Device Name
	Lot Number/Batch Code
	Unique Device Identifier (UDI) as Automatic identification and data capture (AIDC) format
	UDI Device Identifier

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