ENGLISH

SIROS™ Delivery System Technical Description

The information provided below describes the recommended procedural steps for use of the SIROS Delivery System. Read the SIR-Spheres Y-90 resin microspheres Instructions for Use before implanting SIR-Spheres® Y-90 resin microspheres, hereinafter referred to as SIR-Spheres.

The SIROS Delivery System is intended to allow the administration of SIR-Spheres Y-90 resin microspheres.

⚠ CAUTION:

- Patients with diabetes may be at risk of hyperglycemia when 5% dextrose/glucose (D5W) is used to deliver SIR-Spheres and they should be managed using any standard protocols that exist within the context of your treatment site to ensure normoglycemia.

⚠ PRECAUTIONS:

- Time, distance and shielding considerations should be used to minimize exposure to radiation.
- Avoid prolonged stoppages or interruptions to administration process, as this may allow microspheres to settle out of suspension, making them difficult to resuspend and potentially cause blockages.
- Do not over pressurize the system, such as could occur if using syringes smaller than those recommended.

PART NUMBER	DEVICE NAME	DESCRIPTION	PHOTO
SIR-10100	SIROS™ Delivery Dome	Reusable Delivery Dome designed to facilitate the SIR-Spheres administration process and provide beta radiation shielding.	
SIR-10300	SIROS™ Delivery Set	Sterile disposable tubing designed to provide proper routing of dose and contrast media, and connections between the syringes, D-Vial, and the microcatheter set.	
SIR-10200	SIROS™ D-Vial Prep Set	SIROS D-Vial Prep Set is disposable and includes: 1. Sterile D-Vial designed to hold the SIR-Spheres. Supplied with sterile: vent needle (2), filters (2), 21G drawing up needle (1), and blue luer caps (2). 2. D-Vial holder designed to hold the D-Vial. A removeable transport base is provided to support the holder during priming, loading and transport.	

ADDITIONAL ACCESSORIES REQUIRED

- Two 20 mL luer lock syringes filled with non-ionic solution, either 5% dextrose/glucose solution (D5W/G5) or water for injection
- Two 20 mL luer lock syringes filled with non-ionic contrast
- Microcatheter, minimum recommended ID 0.021"

1. SET UP

A. Place the SIROS Delivery Dome on a stable table or cart near the patient. Visually inspect SIROS before use; do not use if damage or wear that could impact function is observed. Ensure that the end marked A is directed toward the microcatheter.

B. Remove the lid from the base by pressing the black unlock button. Set aside lid.


2. TUBING SET UP

A. Remove Delivery Set from sterile pouch.

B. Position Delivery Set on the base with the A line directed toward the microcatheter.

C. Secure Delivery Set lines into their corresponding cleats.

3. PRIMING LINES

- A. Fill two 20 mL syringes with non-ionic solution.
- B. Remove the cap from the D line.
- C. Connect one syringe to the D line and inject non-ionic solution until all air is removed from the D line and fluid exits the A line.

 **CAUTION:** Do not use ionic contrast media with SIR-Spheres.

- D. Fill two 20 mL syringes with non-ionic contrast solution.
- E. Remove the cap from the B line.
- F. Connect one syringe to the B line and inject until all air is removed from the B line and exits the A line.

- G. Remove vented cap from the A line and attach the A line to the microcatheter that is in position in the patient.


 **CAUTION:** Ensure the patient microcatheter is fully connected to the A line before proceeding to step 4.

4. CONNECTING THE D-VIAL

- A. Remove D-Vial holder from its transport base. Place the D-Vial holder into the cupped area on the lower shelf and ensure it is fully seated.

- B. Remove cap from D line on the D-Vial.

- C. Disconnect the D line on the Delivery Set from the red connector and connect it to the D line on the D-Vial.

- D. Remove cap from C line on the D-Vial.
- E. Remove and discard the red connector from the C line on the Delivery Set. Ensure a wet to wet connection is maintained on the C line fittings.
- F. Connect the male luer fitting on the Delivery Set C line to the D-Vial C line.
- G. Place the lid on the base and press down to lock in place. Ensure tubing lines are not pinched.

5. DELIVERY

- A. Deliver the SIR-Spheres suspension in approximately 2 mL aliquots at approximately 5 mL per minute infusion rate using the D Line syringe.

- B. Using the B Line syringe, follow each 2 mL aliquot of SIR-Spheres suspension with a 2 mL aliquot of non-ionic contrast media to facilitate fluoroscopic visualization of flow.

- C. The delivery has been completed when all 20 mL of non-ionic solution has been administered through the SIROS Delivery System. Upon completion, follow disposal procedures below.

6. DISASSEMBLY AND DISPOSAL

- A. Remove syringes from B and D lines. Do not disconnect the patient catheter from the Delivery Set.
- B. Remove the lid from the base by pressing the black unlock button. Set aside lid on the stable table or cart.
- C. Carefully remove catheter from the patient, while still attached to the Delivery Set. Use care when removing and handling the catheter as it may be contaminated with radioactivity. Coil the catheter and wrap within a sterile towel.

- D. Leaving everything attached, lift the D-vial holder out of its cupped area and place it (along with all potentially contaminated procedural accessories) into the waste container for post implant measurements.

 **CAUTION:** Ensure that all items that have come into direct contact with SIR-Spheres are appropriately disposed of as they may be radioactive.

- E. Ensure these items, including 5 mL Luer lock syringe used with Syringe Shield, are disposed of according to facility policy and local regulations.


7. CLEANING OF REUSABLE COMPONENTS

- A. Check SIROS Delivery Dome by using a Geiger Mueller (GM) Detector or wipe test for contamination. Decontaminate if needed.

- B. For cleaning and low level disinfection, use 4 each of a quaternary/alcohol based cleaning/disinfecting wipes such as Super Sani-Cloth. Wipe entire exposed surface and visually inspect to ensure surfaces are clean. Allow surface to remain wet for duration described in the cleaning/disinfecting wipe instructions.

- C. Dry surface with a soft cloth as needed. Do not spray liquid cleaners, use abrasive materials, or sterilize.

SYMBOLS GLOSSARY DEFINITIONS

Symbol	Definition	Symbol	Definition
	Manufacturer		Do not re-sterilize
	Date of Manufacture		Single Use Only. Indicates a medical device that is intended for use on a single patient during a single procedure.
	Consult instructions for use		Product is not made with natural rubber latex
	Caution		Do not use if package is damaged
	Quantity		Keep dry
	Catalog number		Temperature limit
	Lot or batch code		Relative Humidity limit
	Use by date		Ionizing radiation
	Importer		Unique Device Identifier
	Sterile Barrier		Authorized representative in Switzerland
	Protective Barrier		Authorized representative in the European Community
	Non-Sterile		CE mark + Notified Body identification number
	Sterilized using irradiation		Medical Device
	Sterilized using Ethylene Oxide		Caution: Federal law (USA) restricts this device to sale by or on the order of a physician or licensed healthcare practitioner.
	Sterilized using steam		MR Safe
	Serial number		

SIRTeX

Manufacturer
Sirtex Medical Pty Ltd
Shop 6, 207 Pacific Highway
St Leonards, NSW 2065
Australia
Tel +612 9964 8400
Fax +612 9964 8410

Europe, Middle East & Africa
Sirtex Medical Europe GmbH
Joseph-Schumpeter-Allee 33
53227 Bonn
Germany
Tel +49 228 1840 730
Fax +49 228 1840 735

Americas
Sirtex Medical, Inc.
300 Unicorn Park Drive
Woburn, MA 01801
USA
Tel +1 888 4 SIRTEX
Tel +1 888 4 747839
Fax +1 781 721 3880

Asia Pacific
Sirtex Medical Singapore Pte Ltd
50 Science Park Road
#01-01 The Kendall
Singapore Science Park II
Singapore 117406
Tel +65 6800 6500
Fax +65 6778 9424

EC REP

Donawa Lifescience
Consulting Srl
Piazza Albania, 10
00153 Rome, Italy

CH REP

MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

Sirtex Germany Manufacturing GmbH
Branch Frankfurt am Main
Industriepark Hoechst, building G808
65926 Frankfurt am Main
Germany

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