ENGLISH

# SIROS™ Dose Preparation Technical Description

The information and steps described below are recommended for the preparation of the SIR-Spheres® Y-90 resin microspheres within a Nuclear Medicine department. Read the SIR-Spheres Y-90 resin microspheres Instructions for Use and Training Manual before preparing SIR-Spheres Y-90 resin microspheres, hereinafter referred to as SIR-Spheres. The SIROS Delivery System is intended to allow the administration of SIR-Spheres Y-90 resin microspheres.

## ⚠ PRECAUTIONS

- Time, distance and shielding considerations should be used to minimize exposure to radiation.
- Ensure that the distance between any puncture holes in the septum of the D-Vial is at least 2 mm apart.

PART NUMBER	DEVICE NAME	DESCRIPTION	PHOTO
SIR-Y001	SIR-Spheres® Y-90 resin microspheres	SIR-Spheres Y-90 resin microspheres glass vial in Lead Pot	
SIR-10200	SIROS™ D-Vial Prep-Set	SIROS D-Vial Prep Set Consists of: <ul style="list-style-type: none"> <li>• D-Vial transport base</li> <li>• D-Vial holder</li> <li>• D-Vial</li> <li>• 21 G drawing up needle</li> <li>• Two 25 G vent needles</li> <li>• Two 0.22 micron needle filters</li> <li>• Two blue caps</li> </ul>	
SIR-S001	Syringe Shield	Acrylic Syringe Shield	

## ADDITIONAL ACCESSORIES REQUIRED

- Forceps
- Alcohol swabs
- One 20 mL luer lock syringe
- One 5 mL luer lock syringe
- 20 mL of 5% dextrose / glucose solution (D5W/G5) or water for injection

## 1. UNPACKING PROCESS

**A.** Unpack SIR-Spheres lead pot from the plastic shipping pail and place it on the bench top.

**B.** Insert the D-Vial holder into the transport base and place it on the prep surface area.

**C.** Remove D-Vial from pouch and place in D-Vial holder.

**D.** Remove the aluminum seal from the center of the D-Vial with forceps and wipe the rubber septum with an alcohol swab.

## 2. PRIMING PROCESS

**A.** Connect the 0.22 micron filter to the 25 G vent needle. Insert it through the rubber septum of the D-Vial to create a vent.

**B.** Ensure vent needle tip is above the fill level mark.

**C.** Remove the blue cap from D line.

**D.** Attach a syringe filled with at least 10 mL of D5W or sterile water for injection to prime the D line. Fill D-Vial to fill level mark. Ensure there is no air in the D line.

**E.** Disconnect priming syringe from D line; attach a new blue cap to the D line connector.


**F.** Remove blue cap from C line.

**G.** Connect same priming syringe to the C line.

**H.** Slowly draw fluid into the syringe from the vial to prime the C line until fluid is at the prime level mark. Ensure there is no air in the C line. **Do not draw below the prime level mark.**

**I.** Disconnect syringe from the C line; attach a new blue cap to the C line connector.

### 3. DOSE DRAW PROCESS

**A.** Invert the lead pot and shake vigorously before opening to re-suspend the SIR-Spheres, which will have settled during shipping.

**B.** Quickly open the lead pot and remove the glass shipping vial using forceps.

**C.** Determine the total activity of SIR-Spheres in the glass shipping vial using an appropriate ion chamber (dose calibrator) and then return the glass shipping vial to the lead pot.

**D.** Determine the volume of SIR-Spheres suspension that needs to be withdrawn from the glass shipping vial to provide the intended patient-specific activity of SIR-Spheres.

**E.** Partially peel back the aluminum seal of the SIR-Spheres glass shipping vial and wipe the septum with an alcohol swab.

**F.** Connect the 0.22 micron filter to the 25 G vent needle. Insert it through the rubber septum of the shipping vial to create a vent. Ensure that the tip of the needle is well clear of the contents.

**G.** Attach the 21 G drawing up needle to a 5 mL Luer lock syringe and place it in the acrylic syringe shield.

**H.** Using the shielded syringe and 21 G needle, puncture the septum of the glass shipping vial and quickly draw back and forth at least six times in order to re-suspend the SIR-Spheres thoroughly.

**I.** Quickly withdraw the determined volume of SIR-Spheres suspension that will provide the intended patient specific activity. Prior to removing syringe needle from Shipping Vial, draw some air through needle into syringe to draw spheres within needle up into shielded syringe.

**J.** Withdraw the needle from the septum and re-cap the needle using forceps. Set aside on the bench top.

**K.** Using forceps, swirl the glass shipping vial to re-suspend the microspheres and measure the activity remaining in the shipping vial with the dose calibrator.

**L.** Replace glass shipping vial into the lead pot.

**M.** Subtract the activity remaining in the shipping vial from the starting total activity, to determine the amount of activity that has been drawn up into the 5 mL syringe.

**N.** If the amount of activity that has been drawn up into the 5 mL syringe is not correct, transfer the SIR-Spheres suspension back into the glass shipping vial and repeat the steps above to obtain the prescribed level of activity.

**O.** Once the correct activity has been obtained, remove the vent needle from the shipping vial, and secure the lead pot cover.

**P.** Transfer the SIR-Spheres from the 5 mL syringe into the vented D-Vial.

**Q.** If the volume of fluid in the D-Vial does not reach the fill level mark, use a 25 G needle to add D5W or sterile water for injection until the fluid reaches the fill level mark. The volume between the prime level mark and fill level mark is 5 mL.

**Do not exceed the fill level mark.**


**R.** Remove all needles from the D-Vial septum.

**S.** Screw on the D-Vial Holder Cap.

**T.** The patient-specific activity of SIR-Spheres is now ready for transport to the angiography suite in which the implantation will be performed.

## SYMBOLS GLOSSARY DEFINITIONS

Symbol	Definition	Symbol	Definition
	Manufacturer		Do not re-sterilize
	Date of Manufacture		Single Use Only. Indicates a medical device that is intended for use on a single patient during a single procedure.
	Consult instructions for use		Product is not made with natural rubber latex
	Caution		Do not use if package is damaged
	Quantity		Keep dry
	Catalog number		Temperature limit
	Lot or batch code		Relative Humidity limit
	Use by date		Ionizing radiation
	Importer		Unique Device Identifier
	Sterile Barrier		Authorized representative in Switzerland
	Protective Barrier		Authorized representative in the European Community
	Non-Sterile		CE mark + Notified Body identification number
	Sterilized using irradiation		Medical Device
	Sterilized using Ethylene Oxide		Caution: Federal law (USA) restricts this device to sale by or on the order of a physician or licensed healthcare practitioner.
	Sterilized using steam		MR Safe
	Serial number		

**Manufacturer**  
Sirtex Medical Pty Ltd  
Shop 6, 207 Pacific Highway  
St Leonards, NSW 2065  
Australia  
Tel +612 9964 8400  
Fax +612 9964 8410

**Europe, Middle East & Africa**  
Sirtex Medical Europe GmbH  
Joseph-Schumpeter-Allee 33  
53227 Bonn  
Germany  
Tel +49 228 1840 730  
Fax +49 228 1840 735

**Americas**  
Sirtex Medical, Inc.  
300 Unicorn Park Drive  
Woburn, MA 01801  
USA  
Tel +1 888 4 SIRTEX  
Tel +1 888 4 747839  
Fax +1 781 721 3880

**Asia Pacific**  
Sirtex Medical Singapore Pte Ltd  
50 Science Park Road  
#01-01 The Kendall  
Singapore Science Park II  
Singapore 117406  
Tel +65 6800 6500  
Fax +65 6778 9424

Donawa Lifescience  
Consulting Srl  
Piazza Albania, 10  
00153 Rome, Italy

MedEnvoy Switzerland  
Gotthardstrasse 28  
6302 Zug  
Switzerland

Sirtex Germany Manufacturing GmbH  
Branch Frankfurt am Main  
Industriepark Hoechst, building G808  
65926 Frankfurt am Main  
Germany

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