

# SIROS<sup>™</sup>

## Control your Y-90 delivery with confidence

### Driving Innovation for Y-90 Delivery

SIROS is designed to be simple to use yet versatile enough to allow expanded options for patient-tailored delivery



#### Visual

- Visually controlled administration enables precision



#### Intuitive

- Intuitive design for easy set-up, fluid administration and quick disposal



#### Versatile


- Enhanced control allows for more options to tailor sphere delivery



To see a SIROS demonstration, scan here.

# SIROS™

## Intuitive Y-90 Delivery



### Locking cover

Shields and safely contains the Y-90 during delivery.


### Needleless vial

The D-Vial is designed to eliminate touch points and allow more controlled delivery options that may result in a more even distribution of spheres exiting the vial.



### Delivery Set

The peel and place tubing set is designed for quick and simple set-up and disposal.



## Simple Set-up. As easy as 1-2-3

- 1 Position**  
tubing set on base
- 2 Place**  
needleless D-Vial holder
- 3 Snap**  
on cover to lock and secure SIR-Spheres



For more information,  
scan here.

Part Number	Part Name
SIR-10100	SIROS Delivery Dome
SIR-10300	SIROS Delivery Set
SIR-10200	SIROS D-Vial Prep Set

**Caution:** Federal (USA) law restricts this device to sale by or on the order of a physician. Indications for Use: SIR-Spheres Y-90 resin microspheres are indicated for the treatment of unresectable metastatic liver tumors from primary colorectal cancer with adjuvant intrahepatic artery chemotherapy (IHAC) of FUDR (Floxuridine). **Consult the Instructions for Use (www.sirtex.com) for a complete listing of indications, contraindications, side effects, warnings, and precautions.**

### Manufacturer

Sirtex Medical Pty Ltd Shop 6, 207 Pacific Highway St Leonards, NSW 2065 Australia Tel: +61 2 9964 8400 Fax: +61 2 9964 8410  
Americas Sirtex Medical Inc. 300 Unicorn Park Drive Woburn, MA 01801 USA Tel: +1 888 474 7839

www.sirtex.com @SirtexMedical company/sirtex-medical-limited @SirtexMedicalUS

SIR-Spheres® is a registered trademark of Sirtex SIR-Spheres Pty Ltd. ©2023 Sirtex Medical Inc. APM-US-009-05-23

