

FLEXdose Delivery Programme

The Flexibility of Delivery Options


Ordering deadlines 16:00 CET

◆ Wednesday of preceding week

◆◆ Friday of preceding week



Delivery Day



Check availability of delivery options in your country

Delivery dependent on license

* 3 GBq +/- 10% Dose calibrated at 23:00 GMT. Expires 24 hours after calibration day and time

All activity based on 07:00 GMT

FLEXdose Delivery Programme Options

Order Deadline (16:00 CET preceding week to delivery W-1)	Delivery Day	Delivered Activity <i>(estimated activity based on 07:00 GMT)</i>	Days Pre-Calibration	Calibration Day*
	Monday	6.0 GBq	2	Wednesday
		4.6 GBq	1	Tuesday
Wednesday W-1	Tuesday	4.6 GBq	1	Wednesday
		3.6 GBq	0	Tuesday
	Wednesday	3.6 GBq	0	Wednesday
Friday W-1	Wednesday	10.0 GBq	4	Sunday
		4.6 GBq	1	Thursday
	Thursday	7.8 GBq	3	Sunday
		3.6 GBq	0	Thursday
	Friday	6.0 GBq	2	Sunday

* 3 GBq +/- 10% at 23:00 GMT. Expires 24 hours after calibration day and time



Scan here for
Online Order Form

<https://sirtex.com/eu/products/product-ordering/online-order-form-eu/>

Sirtex Medical Europe GmbH
Joseph-Schumpeter-Allee 33, 53227 Bonn, Germany

SIR-Spheres® is a registered trademark of Sirtex SIR-Spheres Pty Ltd

©2022 Sirtex Medical Europe GmbH APM-EMEA-015-07-20 V3