

SIROS™

Control your Y-90 delivery with confidence

Driving Innovation for Y-90 Delivery

SIROS™ is easy to use, visual and adaptable for patient-tailored delivery



Visual

- Full visualization from start to finish



Intuitive

- Intuitive design for easy set-up and streamlined administration



Adaptable

- Adapt your injection for controlled patient-tailored treatment



To see a SIROS demonstration, scan here.



SIROS™

Intuitive Y-90 Delivery



Lockable cover

Shields and safely contains the Y-90 during delivery.



D-Vial

Needleless vial to eliminate touch points and administer safely with even distribution



Delivery Set

The ready-to-use tubing set is designed for streamlined set-up and disposal.



Intuitive Set-up. As easy as 1-2-3

- 1 Position**
tubing set on base
- 2 Place**
needleless D-Vial holder
- 3 Snap**
on cover to lock and secure SIR-Spheres



For more information,
scan here.

Part Number	Part Name
SIR-10100	SIROS Delivery Dome
SIR-10300	SIROS Delivery Set
SIR-10200	SIROS D-Vial Prep Set

Intended Purpose: SIR-Spheres Y-90 resin microspheres are intended for Brachytherapy and are permanently implanted via the hepatic artery to deliver therapy to hepatic tumors. **Indications For Use:** SIR-Spheres Y-90 resin microspheres are indicated for the treatment of: unresectable hepatocellular carcinoma (HCC) or unresectable metastatic liver tumors from primary colorectal cancer (mCRC) in patients refractory to or intolerant of chemotherapy. **Consult the Instructions for Use for a complete listing of indications, contraindications, side effects, warnings, and precautions.** September 2024

Manufacturer

Sirtex Medical Pty Ltd Shop 6, 207 Pacific Highway St Leonards, NSW 2065 Australia Tel: +61 2 9964 8400 Fax: +61 2 9964 8410
Europe, Middle East & Africa Sirtex Medical Europe GmbH Joseph-Schumpeter-Allee 33 53227 Bonn Germany Tel: +49 228 1840 730

www.sirtex.com @SirtexMedical company/sirtex-medical-limited

SIR-Spheres® is a registered trademark of Sirtex SIR-Spheres Pty Ltd. ©2024 Sirtex Medical Europe GmbH APM-EMEA-006-11-23 V4

SIRTeX